

## Adjustable Stroller Handle - Aspire Rehab RX Wheelchairs



- Ergonomic hand grip for carers - especially helpful for tilt and recline movements.
- Height and width adjustable to suit all Rehab RX/RX JNR chair sizes.
- Padded push handle with brakes.
- Pivots up and down via button press positioning.
- The Stroller Handle together with brakes, can be pivoted up and down to suit the carers positioning needs.
- The Adjustable Stroller Handle is mounted on the REHAB RX/RX JNR backrest, replacing the standard handle canes. The attendant brakes are removed from the standard handle canes and re-attached to the Adjustable Stroller Handle (mounting instructions are included in the product carton).



### Warranty Information

- 12 Months

| Code | Display Name | RRP | Box Quantity |
|------|--------------|-----|--------------|
|------|--------------|-----|--------------|

WCA848590 Adjustable Stroller Handle - Aspire Rehab RX Wheelchairs 334 0