

REHAB RX TILT-IN-SPACE WHEELCHAIR

Highly configurable and geared towards the rehabilitation environment, inclusive of pressure care cushioning surfaces, the Rehab RX is designed for users with more complex needs.









Brake Options for Rear and Carer: Dual systems, including push-to-lock and drum brakes, provide users with independence and carers with added safety control

Customisable Seating Dimensions: Adjustable seat width, depth, backrest height, armrest position, and headrest placement allow for tailored postural support to meet individual user needs

Removable Rear Wheels and Height-adjustable Castors: Enhance manoeuvrability, simplify storage in tight spaces, and support vehicle loading

Elevating and Swing-away Legrests: Assist with leg positioning to improve blood circulation, support foot and ankle posture, and enable easier side or front transfers

Anti-tippers Included as Standard: Increase rear stability when tilting or navigating slopes and ramps

ACCESSORIES

IV Pole WCA848280 Safety Lap Belt WCA848300Tray Table WCA848570 Stump Support Various sizes



ADJUSTMENTS



Backrest Height



Seat and Backrest Depth



Elevating Legrest



Tilt-in-space 26° Tilt | 30° Recline

REHAB RX WHEELCHAIR	Seat Width	Seat Depth	Front Seat Height	Rear Seat Height	Overall Width	Backrest Height	Armrest Height	Min Chair Weight	SWL
	mm / in	mm / in	mm / in	mm / in	mm / in	mm / in	mm / in	kg / lbs	kg / st
MWS449710	410 - 460/	440 - 530/	410 - 440/	400 - 475/	690/	460 - 710/	310 - 460/	41/	180/
	16 - 18	17.3 - 20.8	16.1 - 17.3	15.8 - 18.7	27.2	18 - 28	12.2 - 18.1	90.4	28.3
MWS449720	460 - 510/	440 - 530/	410 - 440/	400 - 475/	740/	460 - 710/	310 - 460/	43/	180/
	18 - 20	17.3 - 20.8	16.1 - 17.3	15.8 - 18.7	29.1	18 - 28	12.2 - 18.1	94.8	28.3
MWS449740	510 - 560/	470 - 560/	410 - 440/	400 - 475/	790/	460 - 710/	310 - 460/	47.2/	180/
	20 - 22	18.5 - 22.1	16.1 - 17.3	15.8 - 18.7	31.1	18 - 28	12.2 - 18.1	104.1	28.3

Note: Range of motion dimensions are accurate to a +/-10mm tolerance