

## Mortuary Spreader Bar 363kg



- To be used in conjunction with the Morgue Sling (4 bands) to allow optimal transfers in morgue applications
- The carry bar has a SWL of 363kg
- Accessory is suitable for use with Fixed Ceiling Lifts
- Compatible with Liko patient lifters & Liko ceiling hoists, Prism ceiling hoists, Aidacare A320, A205 and 200 patient lifters



### Compliance

AS/NZS ISO 10535:2011

### Warranty Information

- 1 Year

sova / dignity in motion

Code	Display Name	Safe Working Load
LIC354197	Mortuary Spreader Bar Lifting strap - Short - 2 eyelets - HDPE Plastic - 50mm wide -	
LIC362300	Mortuary Spreader Bar - 363 kg	363 kg
LIC354198	Mortuary Spreader Bar Lifting strap - Short - 4 eyelets - HDPE Plastic - 50mm wide -	