


Aspire ComfiMotion Fixed Allure Headboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs - improving comfort for users and making life easier for carers.

Fixed headboards attach directly to the bed frame thereby moving with the bed during tilt and height adjustments. **Compatible with ComfiMotion ActiveCare beds only.**

- By moving with the bed this reduces pillow displacement and supports the pillow and head position during tilt
- Pivoting design enhances wall protection by pivoting upon wall contact
- The Bed/Headboard can be positioned in direct proximity to a wall and still perform head-down tilt function
- Built-in handset loops (Double/Queen) keep controls within convenient reach
- Widths to match bed sizes, all 680mm high and 120mm deep

 Contract approved

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Colour
BEA830154	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Long Double	749	1	120 mm	1370 mm	680 mm	-
BEA830156	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Queen	799	1	120 mm	1520 mm	680 mm	Portobello Grey
BEA830152	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - King Single	649	1	120 mm	1070 mm	680 mm	Portobello Grey

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Colour
BEA830150	Aspire ComfiMotion Fixed Headboard - Allure - Portobello Grey - Long Single	599	1	120 mm	920 mm	680 mm	Portobello Grey