

TabTimer MK3 Careousel Automatic Pill Dispenser



TabTimer Careousel Automatic Pill Dispenser The Pivotell® Automatic Pill Dispenser reminds the user by means of an alarm and flashing light when to take their medicines.



- **1-28 Daily Alarm Times** - Accomodates the most complicated of medication regimes
- **29 Compartments** - For easy pharmacy filling (one compartment is left empty)
- **Dual / Multi Dose programming** - Two or more compartments can be programmed to alarm at the same time. For those with many or large pills
- **Start delay function** - Delays activation to allow pharmacists or family members to easily fill in advance away from the user's home
- **Pharmacy label carrier and label case** - supplied to comply with pharmacy labelling requirements
- **Upside down feature** - The medication tray rotation is disabled if the dispenser is left upside down, to avoid pill spillage. The alarm will still sound at a dosage time.
- **Day / time discs supplied** for 1 / 2 / 3 / 4 times daily dosage, and one blank disc for five or more doses per day
- **Battery life** - up to 6-12 months
- **Battery level meter** and warning indicator
- **Secure battery case cover** to deter tampering
- **Four different alarm sounds**
- **Flashing light** for those with hearing impairment
- **Lockable** - only the medication due can be accessed by the user
- **12/24 Hour Time** - The clock can be set to either 12 or 24 hour display.
- **Automatic daylight savings** - The clock is self adjusting for daylight savings, this can be turned on/off as required

- **Security** - An alarm signal (different to the pre set alarm sound) is activated if any fault should arise which effects the rotation of the dispenser
 - **Testing** - The Pivotell dispenser has an inbuilt test function which should be activated periodically
 - **Lock** - The Pivotell dispenser is equipped with a lock to ensure that only the medication which is due is accessible to the user **Weight** - with batteries 480 gms
 - **The Dispenser Box Includes:**
 - 1x Mark 3/11 PivoTell® Automatic Pill Dispenser
 - 2x time discs for 1-4 times a day regime
 - 2x blank time discs for more than 4 times a day
 - Key and battery cover securing screw
 - Prescription booklet and holder to apply pharmacy labels or medication information
 - 4 x AA batteries
 - **Optional Accessories:**
- Filling aids** to assist with accurate filling of the trays.
- **Spare trays, lids and time discs** for off-site filling by family members or pharmacists
 - **Disposable tray system** for off-site filling by pharmacists
 - **A Pivotell Dispenser Tipper** for those with limited strength or wrist dexterity and those who might forget to turn the dispenser upright after dispensing their medication

Code	Display Name
DLG275000	TabTimer MK3 Careousel Automatic Pill Dispenser - Solid White Lid

