

Aspire Clamp-On Bed Support Rail



- Designed with clinical functionality in mind, the Aspire Clamp-On Bed Support Rail enhances patient safety during repositioning, transfers, and transitions to and from mobility aids like wheelchairs
- The clamp design ensures consistent setup throughout different facilities, with simple and seamless training and deployment
- Designed to fit the Aspire Lifestyle Community and AC4 Beds

Warranty Information

• 5 Years



Code	Display Name	Total Width	Total Height

BEA829105 Aspire Clamp-On Bed Support Rail 385 mm 425 mm

