

## Aspire ComfiMotion Freestanding Chesterfield Headboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs - improving comfort for users and making life easier for carers.

Freestanding headboards stand independently behind the bed and can be moved or removed as needed. Contemporary in design, the headboard remains stationary while the bed is raised or lowered, not impacting any wall decoration above the headboard.

- Not attached to the bed - offers full positioning flexibility
- Quick to install or remove without tools
- Softened surface that supports pillows and protects the head
- Enhances room style without limiting bed function
- Supports high-needs equipment integration in care environments
- Available in two modern styles, both at 1140mm high and 110mm deep

✓ Contracted DVA product

✓ Contract approved

### Compliance

IEC 60601-2-52:2009

Code	Display Name	RRP	Box Quantity
BEA830100	Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Long Single	458	1
BEA830110	Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - King Single	471	1
BEA830112	Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Long Double	519	0
BEA830114	Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - Queen	539	0
BEA830116	Aspire ComfiMotion Freestanding Headboard - Chesterfield - Portobello Grey - King	709	1

