

Aspire Raphael X CUSTOM Quattro Lift Recline Chair

Customise the Aspire Raphael X Lift Chair to match your unique style and house
Available in a range of commercial quality fabric for long lasting durability



HEALTHCARE VINYL LUSTRELL RANGER RANGE



- **Soft Touch Fabric** - The soft touch weave and diverse colour palette makes the Manisa range a highly practical and functional solution. These fabrics are comforting, easy to clean and durable, yet still soft to touch.
- **Leather** - A luxurious artisan leather with a soft glove-like texture, will retain its unique suppleness, provided they are regularly cared for. The leather is produced and tested to international standards for properties such as rubbing, flex resistance, light fastness, thickness, and strength.
- **Microsuede** - A Suede like texture combined with high performance characteristics ensures the Eastwood fabric range remains easy to clean and maintain.
- **Healthcare Vinyl** The Lustrell Ranger range is high performance vinyl with HealthGuard Antibacterial and Antimicrobial treatment to assist with infection control and hygiene.
- Custom order lead time varies, approximately 5 - 6 weeks depending on fabric availability.
- Custom orders are subject to terms and conditions, please see page two on the order form below.
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only.
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only.

Raphael X lift recline chair details

Designed for enhanced comfort and customisable support. The Raphael X delivers exceptional comfort with features to suit a diverse range of users. The adjustable backrest includes access for personalised postural adjustment and lateral support,

providing improved upright sitting. Its plush design compliments any modern home environment.

- **Built-in Lateral Supports:** Helps maintain an upright sitting posture, improving comfort and stability
- **Integrated Footrest Extension:** Helps with fatigue and improves lower limb comfort*
- **Adjustable Seat Depth:** Offers 25mm of depth adjustment, allowing taller users to achieve a better-seated position
- **Protective Covers for High-Wear Areas:** Removable covers for the headrest and armrests to keep your chair looking new
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environment
- **All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Classic Backrest Option:** Available as an additional accessory for users who prefer a more traditional support profile *The Petite model is designed without a leg extension, catering to the needs of shorter users.

Warranty Information

- 2 Years

| Code | Display Name | RRP | Box Quantity | Total Width | Total Height | Safe Working Load | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|-------------|--|------|--------------|-------------|--------------|-------------------|--------------|------------|----------------------|------------------|
| CHP227640FA | Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Soft Touch Fabric | 4159 | 1 | 820 mm | 1060 mm | 130 kg | 450 - 475 mm | 470 mm | 450 mm | 450 mm |

| Code | Display Name | RRP | Box Quantity | Total Width | Total Height | Safe Working Load | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|-------------|---|------|-----------------|----------------|-----------------|----------------------|---------------|---------------|-------------------------|---------------------|
| CHP227640LT | Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Leather | 4779 | 1 | 820 mm | 1060 mm | 130 kg | 450 - 475 mm | 470 mm | 450 mm | 450 mm |
| CHP227640MS | Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Microsuede | 4159 | 1 | 820 mm | 1060 mm | 130 kg | 450 - 475 mm | 470 mm | 450 mm | 450 mm |
| CHP227640VL | Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Healthcare Vinyl | 4159 | 1 | 820 mm | 1060 mm | 130 kg | 450 - 475 mm | 470 mm | 450 mm | 450 mm |
| CHP227650FA | Raphael X CUSTOM Quattro Lift Recline Chair - SMALL Soft Touch Fabric | 4159 | 1 | 850 mm | 1060 mm | 130 kg | 480 - 505 mm | 500 mm | 450 mm | 450 mm |
| CHP227650LT | Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Leather | 4779 | 1 | 850 mm | 1060 mm | 130 kg | 480 - 505 mm | 500 mm | 450 mm | 450 mm |
| CHP227650MS | Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Microsuede | 4159 | 1 | 850 mm | 1060 mm | 130 kg | 480 - 505 mm | 500 mm | 450 mm | 450 mm |
| CHP227650VL | Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Healthcare Vinyl | 4159 | 1 | 850 mm | 1060 mm | 130 kg | 480 - 505 mm | 500 mm | 450 mm | 450 mm |
| CHP227660FA | Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM Soft Touch Fabric | 4159 | 1 | 875 mm | 1100 mm | 130 kg | 500 - 525 mm | 525 mm | 470 mm | 470 mm |
| CHP227660LT | Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Leather | 4779 | 1 | 875 mm | 1100 mm | 130 kg | 500 - 525 mm | 525 mm | 470 mm | 470 mm |
| CHP227660MS | Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Microsuede | 4159 | 1 | 875 mm | 1100 mm | 130 kg | 500 - 525 mm | 525 mm | 470 mm | 470 mm |
| CHP227660VL | Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Healthcare Vinyl | 4159 | 1 | 875 mm | 1100 mm | 130 kg | 500 - 525 mm | 525 mm | 470 mm | 470 mm |

| Code | Display Name | RRP | Box Quantity | Total Width | Total Height | Safe Working Load | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|-------------|---|------|-----------------|----------------|-----------------|----------------------|-----------------|---------------|-------------------------|---------------------|
| CHP227670FA | Raphael X CUSTOM Quattro Lift Recline Chair - LARGE Soft Touch Fabric | 4159 | 1 | 900 mm | 1100 mm | 130 kg | 530 - 555 mm | 550 mm | 500 mm | 500 mm |
| CHP227670LT | Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Leather | 4779 | 1 | 900 mm | 1100 mm | 130 kg | 530 - 555 mm | 550 mm | 500 mm | 500 mm |
| CHP227670MS | Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Microsuede | 4159 | 1 | 900 mm | 1100 mm | 130 kg | 530 - 555 mm | 550 mm | 500 mm | 500 mm |
| CHP227670VL | Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Healthcare Vinyl | 4159 | 1 | 900 mm | 1100 mm | 130 kg | 530 - 555 mm | 550 mm | 500 mm | 500 mm |
| CHP227680FA | Raphael X CUSTOM Quattro Lift Recline Chair - MAXI- Soft Touch Fabric | 4369 | 1 | 1020 mm | 1120 mm | 180 kg | 540 - 565 mm | 650 mm | 510 mm | 510 mm |
| CHP227680LT | Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Leather | 4989 | 1 | 1020 mm | 1120 mm | 180 kg | 540 - 565 mm | 650 mm | 510 mm | 510 mm |
| CHP227680MS | Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Microsuede | 4369 | 1 | 1020 mm | 1120 mm | 180 kg | 540 - 565 mm | 650 mm | 510 mm | 510 mm |
| CHP227680VL | Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Healthcare Vinyl | 4369 | 1 | 1020 mm | 1120 mm | 180 kg | 540 - 565 mm | 650 mm | 510 mm | 510 mm |