

Aidacare FL250 Floorline Bed - Telescopic Rails



- Adjustable height rails to assist in the prevention of falls as well as manoeuvring.
- Includes both left/right side rails, and bed height limiting handset.
- Various configurations include 1/3, 2/3 individually or together for full length
- *Bed ends sold separately
- Please note: The minimum bed height changes to 275mm when Telescopic Rails are installed



Code	Display Name	RRP	Box Quantity	Colour
BEB048800	FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Rail - Sublime Teak	1869	1	Sublime Teak
BEB048805	FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Rail - Milano Walnut	1869	1	Milano Walnut
BEB048810	FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Decorative Trim - Sublime Teak	1249	1	Sublime Teak
BEB048815	FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Decorative Trim - Milano Walnut	1249	1	Milano Walnut
BEB048820	FL250 Telescopic Rails - 2x2/3 Rail + 2x1/3 Decorative Trim - Sublime Teak	1349	1	Sublime Teak
BEB048825	FL250 Telescopic Rails - 2x2/3 Rail + 2x1/3 Decorative Trim - Milano Walnut	1349	1	Milano Walnut
BEB048830	Bed End Set - FL250 - Single - ComoR (Sublime Teak) - For Telescopic Rails	-	-	Sublime Teak

Code	Display Name	RRP	Box Quantity	Colour
BEB048835	Bed End Set - FL250 - Single - ComoR (Milano Walnut) - For Telescopic Rails	-	-	Milano Walnut
BEB048840	Bed End Set - FL250 - King Single - ComoR (Sublime Teak) - For Telescopic Rails	-	-	Sublime Teak
BEB048845	Bed End Set - FL250 - King Single - Como R (Milano Walnut) - For Telescopic Rails	-	-	Milano Walnut