

Aidacare FLX Floorline Bed



- Ultra Low Height: Market leading 100mm low height whilst maintaining full ground clearance
- Pull Apart Capability: The bed can be disassembled into two sections, offering improved storage and portability.
- Attendant Control Panel: Allows caregivers functional control and the ability to lock out certain bed functions from being operated via the handset
- Underbed Clearance: Full clearance under the bed minimises the risk of cord impingement and protects flooring from damage
- Custom Height Function: Set the perfect height for the individual user. Caters to a large range of user heights
- SmartDescent: Bed lowering speed is reduced once the bed reaches 300mm maximizing carer potential to remove obstacles as the bed lowers to the floor
- Underbed Lighting: Creates a gentle halo around the bed to reduce falls risk and improve carer safety
- Backrest Auto-Regression: Reduces shear forces on the spine and sacral regions. Minimizes abdominal compression when back is raised supporting cardio pulmonary and respiratory function
- One Touch Chair Positioning: Improves socialinteraction, assists respiratory and cardio pulmonary function and creates a comfortable position for sitting, dining and transfers
- Integrated Bed Extension: 100-200mm length extension is easily operated and ensures a large range of user heights will fit comfortably without compromise
- Trendelenburg (and Reverse): Augments patient positoning, comfort and assists with bed mobility*.
- Quick Swap Actuators: Allow for quick and easy inspection and maintenance, minimising downtime
- TiMotion Electronics: Medical grade electronics from a leading global



manufacturer

- Control Box: Easily accessible beneath the back rest panel for efficient servicing
- Quick Connect Bed Ends: Simple and quick to set-up upon bed installation. Easily removed for access or maintenance
- Wall Buffer: Length adjustable to suit room configuration
- Cable Conduit: Allows other cables like air mattress power cords to be safely tucked away, preventing trip hazards and cord damage
- EasyMove Transport System: Facilitates safe no-lift transportation and movement of beds, as well as safe and easy access underneath the bed platform for routine inspection and maintenance
- Width Extension Kit: Ability to increase a single to a king single, ideal for resident changes in a facility
- High Performance Twin Castors x 8: Head end total locking (x2), foot end central locking (x4)

Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset & Cables

| Code | Display Name | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight | Colour |
|-----------|--|----------------|-------------|--------------|-------------------|-------------|---------------|
| BEB048600 | Aidacare FLX Floorline Bed – Single | 2130 - 2330 mm | 964 mm | 100 - 750 mm | 275 kg | 105 kg | Sublime Teak |
| BEB048610 | Aidacare FLX Floorline Bed – King Single | 2130 - 2330 mm | 1124 mm | 100 - 750 mm | 275 kg | 115 kg | Milano Walnut |