

Aidacare Aspire Lifter Leg Protectors



- Adhesive strips designed to wrap around and adhere to the Lifter leg surface, protecting the powder coat against abrasion and scratches caused during normal use
- Each kit contains a surface cleaning wipe and 2 Leg Protector strips – 1 for each leg
- Highly durable and adhesive vinyl protective layer
- Protect against scrapes and chips of the powder coat layer
- Transparent design
- Easily applied and removed
- Covering all Aidacare Aspire Lifters except for the A150F
- Can be applied to other lifter brands - strips can easily be trimmed to the appropriate size, using a blade/scissors



Code	Display Name	RRP	Box Quantity	Total Length	Total Width
GSP390815	Aidacare Aspire Lifter Leg Protectors - A205, A200S & 200S Lifter	88.5	1	650 mm	160 mm
GSP390816	Aidacare Aspire Lifter Leg Protectors - 200 Lifter	88.5	1	780 mm	140 mm
GSP390817	Aidacare Aspire Lifter Leg Protectors - A320 Lifter	88.5	1	780 mm	180 mm