

Oppo 2265 Lumbar Sacral Support - Tapered



- Four removable posterior flexible aluminum stays support the lower back
- Two flexible stays in front help reduce rolling of the support when sitting
- Allows limited flexion and extension
- 25cm high at back, tapered to 17.5cm at front

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	Chest Circumference
FBS301278XL	Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - XL	965 - 1092 mm
FBS301278XXL	Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - XXL	1092 - 1219 mm
FBS301278M	Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - Med	711 - 838 mm
FBS301278	Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - Sml	584 - 711 mm
FBS301278L	Oppo - Sacro Lumbar Suppt w/ Alumn Stays - Tapered - Lge	838 - 965 mm