

Automatic Calendar Clock



- The BQ38 is popular at home as well as in the office
- It features a reliable German movement and an automatic calendar function that self-adjusts for leap years
- This is an ideal clock wherever people are completing forms and need a reminder of the day of the week

Warranty Information

- 1 Year



Support @ Home

Function: Managing time, memory and planning

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard



Code	Display Name
DLG277090	Jadco Analogue Clock with Automatic Flip Calendar