

## Aspire ComfiMotion Fixed Companion Bed



- Matching single bed without any adjustability designed to sit flush when positioned together with the ComfiMotion CARE Bed and ComfiMotion Activ Care Bed
- It facilitates a side-by-side sleeping arrangement at a fixed height of 370 mm, without compromising on the users care and assistance requirements

✓ Contracted DVA product

✓ Contract approved

### Compliance

IEC 60601-2-52:2009

### Warranty Information

- 10 years - Frame
- 1 year - Fabric



Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load
BEB046800	Aspire ComfiMotion Fixed Companion Bed - Portobello Grey - Long Single (Headboard Sep)	1009	1	2030 mm	920 mm	370 mm	180 kg