

Aspire Rembrandt Quattro Lift Recline Chair



The Rembrandt Quattro offers a unique seating experience. With independent control of back, leg, lumbar and headrest functionality, the user has infinite positioning options with the press of a button. Genuine shiatsu massage coupled with optional heating gives users a truly unique and luxurious comfort experience. Zero Gravity tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution.



- **Shiatsu Massage** in backrest provides therapeutic and stress relief for the user
- **Memory Function** allows for a preset favourite position to be set
- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book
- **Integrated Lift and Recline** functions for easy sit-to-stand
- **Extending Legrest** supports more of the lower legs for enhanced comfort
- **Heating Function** within the seat, lumbar and backrest provides warmth and comfort during the colder months
- **Cup Holder** built into the armrest for easy drink storage
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution
- **Wireless Device Charging** built into the armrest for convenience
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day
- **Backlit Hand Control** with integrated USB charging point
- Looking to protect your Electric Lift Recline Chair? Try our [Universal Chair Protector](#)

Warranty Information

- 5 Years - Frame

- 2 Years - Components

| Code | Display Name | RRP | Box Quantity | Total Width | Total Height | Safe Working Load | Colour | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest | Backrest Height | User Height |
|-----------|---|------|--------------|-------------|--------------|-------------------|---------------|------------|------------|----------------------|------------------|-----------------|-------------|
| CHP227550 | Aspire Rembrandt Quattro Lift Recline Chair | 5549 | 1 | 1015 mm | 1095 mm | 159 Kg | Canyon Walnut | 500 mm | 560 mm | 495 mm | 495 mm | 600 mm | 1065 mm |