

## Aidacare Anaesthesia Cart



- The Aidacare Anaesthesia Cart is the flexible solution to your workflow needs
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Combined with its sleek lines, strong construction and easy manoeuvrability, it offers a fully customisable solution to meet individual needs
- Easy roll 125 mm castors
- Packed with standard accessories, it allows for choice of cart height, draw depth, locking option and storage configuration to perfectly compliment your workflow needs
- **Standard Inclusions:** Accessory Tower, Drawer Dividers, Extended Sliding Surface and Waste Bin with Lid
- Choose from an array of accessories to optimise the cart for your facility's treatment needs
- Lock System Options:
  - **Break-away lock with Plastic Seal:** The simple breakaway lock offers easy access, combined with essential security
  - **Digital lock:** A user friendly system where the cart is locked and unlocked using a pre-set code ensuring ease of access and the security required in clinical environments. The digital lock is robust and waterproof
- Standard drawer configuration detailed below
- **Small Anaesthesia Cart:**
  - 1 small 3' drawer (7 cm)
  - 2 medium 6' drawers (16 cm)
  - 1 large 9' drawer (25 cm)
- **Medium Anaesthesia Cart:**



### Accessories

Click on the cart icon to display the range of optional accessories

- 2 small 3' drawers (7 cm)
- 2 medium 6' drawers (16 cm)
- 1 large 9' drawer (25 cm)
- **Large Anaesthesia Cart:**
  - 1 small 3' drawer (7 cm)
  - 3 medium 6' drawers (16 cm)
  - 1 large 9' drawer (25cm)

## Warranty Information

- Structural: 5 Years Warranty
- Digital Lock: 2 Years Warranty
- Accessories: 1 Year Warranty

| Code      | Display Name   | Total Length | Total Width | Total Height | Colour |
|-----------|--|--------------|-------------|--------------|--------|
| WAT909660 | Aidacare Anaesthesia Cart - Small - Break-away Lock  | 790 mm       | 610 mm      | 900 mm       | Green  |
| WAT909630 | Aidacare Anaesthesia Cart - Small - Digital Lock     | 790 mm       | 610 mm      | 900 mm       | Green  |
| WAT909670 | Aidacare Anaesthesia Cart - Medium - Break-away Lock | 790 mm       | 610 mm      | 1000 mm      | Green  |
| WAT909640 | Aidacare Anaesthesia Cart - Medium - Digital Lock    | 790 mm       | 610 mm      | 1000 mm      | Green  |
| WAT909680 | Aidacare Anaesthesia Cart - Large - Break-away Lock  | 790 mm       | 610 mm      | 1100 mm      | Green  |
| WAT909650 | Aidacare Anaesthesia Cart - Large - Digital Lock     | 790 mm       | 610 mm      | 1100 mm      | Green  |