

Dale Abdominal Lumber Support



- Single ply stretch loop elastic (latex-free)
- Two removable aluminium supports
- Supports can be hand-shaped to hip curvature
- Velcro closure anywhere on belt for perfect fit
- Colour: white
- Sizing: available in either 23cm (pictured) or 30cm widths, both in S, M and L sizes

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	Chest Circumference
FBS301272	Dale - Abdominal / Lumbar - Elastic - 23cm wide - fits to 152-191cm - Lge	1520 - 1910 mm
FBS301271	Dale - Abdominal / Lumbar - Elastic - 23cm wide - fits to 114-157cm - Med	1140 - 1570 mm
FBS301270	Dale - Abdominal / Lumbar - Elastic - 23cm wide - fits to 114cm - Sml	1140 mm
FBS301274	Dale - Abdominal / Lumbar - Elastic - 30cm wide - fits to 114-157cm - Med	1140 - 1570 mm
FBS301275	Dale - Abdominal / Lumbar - Elastic - 30cm wide - fits to 152-191cm - Lge	1520 - 1910 mm
FBS301273	Dale - Abdominal / Lumbar - Elastic - 30cm wide - fits to 114cm - Sml	1140 mm