

TabTimer Careousel Dispenser Accessories



- **Spare Tray – MK3:** Includes tray, lid, and time discs. Fits TabTimer Careousel MK3 dispenser.
- **Spare Tray – ADVANCE:** Includes tray, lid, and day/time discs. Fits Careousel ADVANCE dispenser.
- **Dust Cover – MK3:** Protects the MK3 dispenser from dust and damage when not in use.
- **Dust Cover – ADVANCE:** Shields the ADVANCE dispenser for longer-lasting protection.
- **Dosing Ring – MK3:** Set of three double-sided printed discs for accurate dose setup.
- **Dosing Ring – ADVANCE:** Replacement card discs for Careousel ADVANCE dispenser.
- **Spare Key:** Replacement key for TT24-28ADV and TT28-29MK3 models.
- **Shutter – ADVANCE:** Replacement shutter for Careousel ADVANCE dispenser to maintain secure access.
- **Pharmacy Label:** Label kit for identifying patient or medication details on dispensers.
- **Pharmacy Holder:** Secure label holder designed for TabTimer Careousel dispensers.
- **Frosted Lid – MK3:** Replacement frosted lid designed to fit MK3 pill dispensers.
- **Frosted Lid – ADVANCE:** Replacement frosted lid for ADVANCE pill dispensers.
- **Filling Funnel:** Assists with quick and accurate loading of medication trays for MK3/ADVANCE models.
- **Filler Tray – MK3:** Simplifies pharmacy filling process for MK3 dispensers.
- **Filler Tray – ADVANCE:** Designed for easy pharmacy filling of ADVANCE dispensers.



- **Tipper:** Assists users with limited dexterity in dispensing pills safely and easily.
- **Tamper Evident Sticker:** Set of 52 labels that leave a VOID message if tampered with.
- **Battery Cover:** Replacement battery cover for TabTimer Careousel MK3 or ADVANCE models.

Code	Display Name
DLG275080	TabTimer Careousel Pill Dispenser - Tipper
DLG275090	TabTimer - Tamper Evident Sticker
DLG275100	TabTimer MK3 - Spare Tray
DLG275110	TabTimer ADVANCE - Spare Tray
DLG275120	TabTimer Careousel Pill Dispenser - Spare Key
DLG275130	TabTimer ADVANCE - Shutter
DLG275150	TabTimer Careousel Pill Dispenser - Pharmacy Holder
DLG275160	TabTimer MK3 - Frosted Lid
DLG275170	TabTimer ADVANCE - Frosted Lid
DLG275180	TabTimer ADVANCE - Filling Funnel
DLG275190	TabTimer MK3 - Filler Tray
DLG275200	TabTimer ADVANCE - Filler Tray

Code	Display Name
DLG275210	TabTimer MK3 - Dust Cover
DLG275220	TabTimer ADVANCE - Dust Cover
DLG275230	TabTimer MK3 - Dosing Ring
DLG275240	TabTimer ADVANCE - Dosing Ring
DLG275250	TabTimer Careousel Pill Dispenser - Battery Cover