

Kerry General Purpose Sling with Head Support - Mesh



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user
- Made from a lightweight mesh which is ideal for showering or bathing and dries quickly
- Provides excellent back and head support for patients requiring support during transfers
- Ideal for non weight bearing patients
- FIM Scores: 1 and 2



ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS381800	General Purpose Sling with Head Support - Small	250 kg	Mesh
LSS381900	General Purpose Sling with Head Support - Medium	250 kg	Mesh
LSS382000	General Purpose Sling with Head Support - Large	250 kg	Mesh
LSS382100	General Purpose Sling with Head Support - Extra Large	250 kg	Mesh