

NuboMed X3000 Lite Mobile Nursing Cart



- Front facing drawers for conventional operation.
- Hidden and expandable work surface for maximum usage space.
- A variety of accessories and drawer combinations.
- Mechanical lock for security.
- Embedded computer, HD display, design with care.
- Wireless charging (Phone/PDA), RFID user log-in.
- Ideal For: Daily Nursing/Data Recording.



Warranty Information

- 2 Years

NuboMed

Code	Display Name
WAT908180	NuboMed - X3000 - Lite Mobile Nursing Cart