

Bed/Chair Raiser



- Raisers are placed on the floor under the legs of beds or lounge chairs to raise their height
- Ideal for those recovering from hip replacements, knee operations or as a result of ageing

Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_181212202_0103_1_2

Minor Bed – Composite:

05_428_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Height
BEA009700	Bed / Chair Raiser - Blocks - Wooden 75mm (Pkt 4)	212	1	75 mm
BEA009600	Bed/Chair Raiser Cones - Plastics - Pack of 4 - 140mm High	132	1	140 mm