

Diffuser Cushion

- Premium density pressure-relief cushion made from memory foam
- Ideal for wheelchairs or conventional seating



Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

ISO 7176-7:1996

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Pressure Care Rating	Material
GSS340100	Diffuser Cushion	156	1	-	-	-	-
GSS340103	Diffuser Cushion with Waterproof Cover	163	1	420 mm	440 mm	Low Risk	Memory Foam
GSS340104	Diffuser Cushion - Large	177	1	460 mm	490 mm	Low Risk	Memory Foam
GSS340105	Diffuser Cushion with Waterproof Cover - Large	177	1	460 mm	490 mm	Low Risk	Memory Foam
GSS340101	Waterproof Cover - Standard	81	1	420 mm	440 mm	Low Risk	-
GSS340106	Waterproof Cover - Large Size	88	1	460 mm	490 mm	Low Risk	-

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Pressure Care Rating	Material
------	--------------	-----	--------------	--------------	-------------	----------------------	----------