

## Lifting Strap



- To be used in conjunction with the Morgue Sling (4 bands) to allow optimal transfers in morgue applications
- The carry bar has a SWL of 363 kg
- Accessory is suitable for use with Fixed Ceiling Lifts
- Compatible with Liko patient lifters & Liko ceiling hoists, Prism ceiling hoists, Aidacare A320, A205 and 200 patient lifters



### Compliance

AS/NZS ISO 10535:2011

## Warranty Information

- 1 Year

sova / dignity in motion

Code	Display Name	Total Length	Total Width	Total Height
LIC354199	Mortuary Spreader Bar with eyelets	2000 mm	80 mm	3 mm
LIC362400	Lifting strap - Short - 2 eyelets - HDPE Plastic - 50mm wide	1250 mm	50 mm	-
LIC362410	Lifting strap - Long - 2 eyelets - HDPE Plastic - 50mm wide	1475 mm	50 mm	-