

## Lifting Strap



- To be used in conjunction with the Morgue Sling (4 bands) to allow optimal transfers in morgue applications
- The carry bar has a SWL of 363 kg
- Accessory is suitable for use with Fixed Ceiling Lifts
- Compatible with Liko patient lifters & Liko ceiling hoists, Prism ceiling hoists, Aidacare A320, A205 and 200 patient lifters



### Compliance

AS/NZS ISO 10535:2011

## Warranty Information

- 1 Year

sova / dignity in motion

| Code      | Display Name   | Total Length | Total Width | Total Height |
|-----------|--|--------------|-------------|--------------|
| LIC354199 | Mortuary Spreader Bar with eyelets                           | 2000 mm      | 80 mm       | 3 mm         |
| LIC362400 | Lifting strap - Short - 2 eyelets - HDPE Plastic - 50mm wide | 1250 mm      | 50 mm       | -            |
| LIC362410 | Lifting strap - Long - 2 eyelets - HDPE Plastic - 50mm wide  | 1475 mm      | 50 mm       | -            |