

Aspire ComfiMotion Allure Footboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs—improving comfort for users and making life easier for carers.



The footboard completes the bed in both style and function, compatible with **ComfiMotion ActiveCare beds only.**

- Helps complete the bed end and keep bedding in place during profiling, reducing the need for carers to continually re-adjust bedding
- Provides a practical accessory mounting surface (e.g. mattress pump)
- Portobello grey upholstery to match bed and headboards
- Secure and sturdy bolt-on attachment to bed frame
- Widths to match bed sizes, all 610mm high and 60mm deep

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Colour
BEA830166	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Queen	400	1	60 mm	1520 mm	610 mm	Portobello Grey
BEA830164	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Long Double	375	1	60 mm	1370 mm	610 mm	Portobello Grey
BEA830162	Aspire ComfiMotion Footboard - Allure - Portobello Grey - King Single	350	1	60 mm	1070 mm	610 mm	Portobello Grey
BEA830160	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Long Single	325	1	60 mm	920 mm	610 mm	Portobello Grey

