

Aspire ComfiMotion Memory Foam Mattress - Pressure Care



- The ComfiMotion Memory Foam mattress contours to the user's body, providing personalised comfort and transitional pressure relief
- Breathable, 4-Way Stretch Waterproof Cover to reduce shear, and friction and maintain skin microclimate
- Medical Grade Waterfall Cover is fully welded over the zipper to maximise infection control
- Triple-layered foam:
 - Top Layer: Immersive visco-elastic
 - Middle Layer: High resilience foam for support and pressure redistribution
 - Base Layer: Strong base support for maximum durability
- Integrated side walls increase roll resistance at the bed edge and assist users, partners, or carers with extra rigidity to ease bed entry and exit
- Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief
- Superior breathability: Up to 8 x more breathable than conventional memory foams
- Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria
- Eco-friendly GECA Foam is produced using cleaner chemicals and production technology

- ✔ Contracted DVA product
- ✔ Contract approved

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load
BEM469930	Aspire ComfiMotion Memory Foam Mattress Pressure Care - Long Single	1149	1	2030 mm	920 mm	200 mm	180 kg

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load
BEM469940	Aspire ComfiMotion Memory Foam Mattress Pressure Care - King Single	1229	1	2030 mm	1070 mm	200 mm	180 kg
BEM469942	Aspire ComfiMotion Memory Foam Mattress Pressure Care - Long Double	1539	0	2030 mm	1370 mm	200 mm	360 kg
BEM469944	Aspire ComfiMotion Memory Foam Mattress Pressure - Queen	1659	0	2030 mm	1520 mm	200 mm	360 kg