

Kineris Air Cushion



Accurate Pressure Adjustment: Manometer enables exact pressure adjustment every time. Easy for carers.



- **Single-valve:** For simple inflation and maintenance.
- **Reduced Friction and Shearing:** Cell mobility minimises friction and shearing forces.
- **Optimal Height:** Available in low and high profiles.
- **Comfortable and Breathable:** The removable cover is stretchable, and breathable, and aids body temperature regulation during prolonged sitting.
- **Stays in Place:** The inner cover features a non-slip base to keep the cushion secure.
- **Versatile Sizing:** Available in various sizes to accommodate different needs.

Warranty Information

- 3 years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC627100	Kineris Air Cushion - Low Profile - 16.5 x 16.5in	16.5 in	16.5 in	2.75 in	150 kg
PTC627110	Kineris Air Cushion - Low Profile - 16.5 x 18.5in	16.5 in	18.5 in	2.75 in	150 kg
PTC627120	Kineris Air Cushion - Low Profile - 18.5 x 18.5in	18.5 in	18.5 in	2.75 in	150 kg
PTC627200	Kineris Air Cushion - High Profile - 16.5 x 16.5in	16.5 in	16.5 in	4 in	150 kg
PTC627210	Kineris Air Cushion - High Profile - 16.5 x 18.5in	16.5 in	18.5 in	4 in	150 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC627220	Kineris Air Cushion - High Profile - 18.5 x 18.5in	18.5 in	18.5 in	4 in	150 kg