

JAY - Care Backrest



- Designed to accommodate a client who has a fixed kyphotic, the JAY-CARE backrest provides solid back support to accommodate and slow-down back deformity
- The contoured design and 39° angle adjustability help further accommodate the multiple needs of the kyphotic client and make the JAY-CARE suitable to a client who has changing postural needs over time
- Provides mild to moderate lateral support through the built-in lateral contour of the shell
- Optional foam wedges provide for a more customised fit
- All models include: Jay AirExchange cover and universal mounting hardware

Warranty Information

- 2 Years Backrest
- 6 Months Cover



ndis

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_417_0105_1_2

05_419_0105_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

| Code | Display Name | Seat Width | Backrest Width | Backrest Height | Safe Working Load |
|-----------|-----------------|------------|----------------|-----------------|-------------------|
| PTB532530 | JAY Care - 16in | 15 - 16 in | 16 in | 16.5 in | 113 kg |
| PTB532520 | JAY Care - 18in | 17 - 18 in | 18 in | 17.25 in | 113 kg |
| PTB532515 | JAY Care - 20in | 19 - 20 in | 20 in | 19 in | 113 kg |