


Aspire Hybrid Comfort Mattress



The HYBRID COMFORT is an Acute mattress system offering foam-in-cell technology enabling flexible operation as a high comfort foam static support surface or a powered-up alternating pressure support surface. Specifically re-engineered for increased patient comfort, foam layering now effectively supports patient weight ranging from 30 - 230 kg.

- **Immersive Visco-elastic Memory Foam Surface:** Top layer is now twice as thick to offer significant improvement in user comfort
- **Intelligent Hose Management:** When not in use the air-hose can be stored internally, reducing infection control and OH&S risks
- **Integrated Shear Reduction:** Memory foam layer moves independently of the cell surface reducing shear when bed is adjusted
- **Optimised Heel Protection:** Tapered heel zone redistributes pressure to lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear
- **Innovative pump technology:** Reduces noise and vibration for uninterrupted sleep

- ✓ Contracted DVA product
- ✓ Contract approved

Prescribing Therapists: Click here to download our [Clinical Assistance Guide](#)

Warranty Information

- 5 Years Foam
- 4 Years Cover
- 2 Years Pump/cells

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load
PTM638005	Aspire Hybrid Comfort Mattress - Single	2149	1	1980mm	895mm	170mm	30-230 kg
PTM638006	Aspire Hybrid Comfort Mattress - King Single	2399	1	1980mm	1050mm	170mm	30-230 kg
PTM638030	Digital Pump	759	1	320 mm	230 mm	130 mm	-
PTM638007	Aspire Hybrid Comfort Mattress - Bariatric	3099	1	1980 mm	1050 mm	170 mm	30-380 kg
PTM638032	Alternating Pump - Bariatric	739	1	-	-	-	-