

Aspire Go Turner



The Aspire Go Turner is a proven sit-to-sit transfer aid for quick and easy manual handling. This compact folding patient mover has a very narrow footprint, making it the perfect solution for transferring patients in tight spaces. The pull apart and folding design reduces storage space requirements and helps to transport the product if needed.



- Ergonomic handlebar
- Sling retention clips - Large range of sling belts (S to XL) for improved support and positioning
- Angle and height adjustable knee pads offer comfort during transport
- Central locking footbrake
- Spacious footplate
- 6 x wheels for improved stability

| Code | Display Name | RRP | Box Quantity | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight | Waist/Hip Circumference |
|-----------|--|------|--------------|--------------|-------------|--------------|-------------------|-------------|-------------------------|
| LSS390585 | Aspire Go Turner - Sit-to-Sit Transfer Aid | 2309 | 1 | 685 mm | 570 mm | 1205 mm | 150 kg | 18.9 kg | - |
| LSS391132 | Aspire Go Turner Standing Belt - Small | 198 | 1 | - | - | - | 150 kg | - | 780 - 980 mm |
| LSS391133 | Aspire Go Turner Standing Belt - Medium | 198 | 1 | - | - | - | 150 kg | - | 930 - 1180 mm |
| LSS391134 | Aspire Go Turner Standing Belt - Large | 198 | 1 | - | - | - | 150 kg | - | 1130 - 1380 mm |
| LSS391135 | Aspire Go Turner Standing Belt - Xlarge | 198 | 1 | - | - | - | 150 kg | - | 1330 - 1530 mm |

| Code | Display Name | RRP | Box Quantity | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight | Waist/Hip Circumference |
|------|--------------|-----|--------------|--------------|-------------|--------------|-------------------|-------------|-------------------------|
|------|--------------|-----|--------------|--------------|-------------|--------------|-------------------|-------------|-------------------------|