

## **SitnStand Lift Assist - Classic**



Upgrade your home and existing chair with our Lift Assist Classic, designed to enhance mobility and independence. This versatile device offers reliable assistance for individuals who require support when transitioning between sitting and standing positions at a fraction of the cost of a lift chair.

- Sturdy Construction: Built with durable materials for long-lasting performance and stability.
- Simple Installation: Easily install without the need for complicated tools or professional assistance.
- Watch the informational video here.
- Smooth Operation: Experience smooth and effortless transitions from sitting to standing with the push of a button.
- Comfortable Design: The lift assist is designed with user comfort in mind.
- Compact Size: Its compact size allows for easy placement in various areas of your home, maximising convenience.
- Weight Allowance: Up to 266lbs/121kg.
- Rechargeable Battery: Last for 4-7 days.
- Washable Seat Cover: Easy to remove.

## Warranty Information

• 2 years

Code Display Name Safe Working Load

LSS403710 SitnStand Lift Assist - Classic 121 kg



