

BC2 Waterfall Day Chair



- High back orthopaedic ergonomic chair
- Three segment padded back cushion, padded arm rests and seat for added comfort
- All cushions are removable from welded one piece shell for easy cleaning
- Fully adjustable seat height
- Powder coated steel frame
- Range of upholstered finishes and colours available
- Please specify colour at time of order

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_188800051_0103_1_2

Click [here](#) for more information

Warranty Information

- 10 Years Structural
- 1 Year Fabric

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Total Height	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP186050	Chair Tray and Attachments	381	-	-	-	-	-	-	-	-
CHP196598	BC2 Waterfall Day Chair - Vinyl - Kiwi	899	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Vinyl_Kiwi
CHP196604	BC2 Waterfall Day Chair - Vinyl - Ink	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Vinyl Ink
CHP196606	BC2 Waterfall Day Chair - Vinyl - Latte	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Vinyl Latte

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Total Height	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP196608	BC2 Waterfall Day Chair - Fabric - Tweed Charcoal	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Tweed Charcoal
CHP196610	BC2 Waterfall Day Chair - Fabric - Tweed Oatmeal	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Tweed Oatmeal
CHP196612	BC2 Waterfall Day Chair - Fabric - Tweed Denim	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Tweed Denim
CHP196614	BC2 Waterfall Day Chair - Fabric - Tweed Stone	919	1	490 mm	510 mm	1020 -1190 mm	680 mm	400 - 550 mm	160 kg	Tweed Stone