

BC1 High Back Day Chair



- Stylish orthopaedic ergonomic chair
- Padded armrests and padded fully welded one piece shell
- Fully adjustable seat height
- Powder coated steel frame
- Range of upholstered finishes and colours available
- Please specify colour at time of order

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_188800051_0103_1_2

Click [here](#) for more information

Warranty Information

- 10 Years Structural
- 1 Year Fabric

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Colour
CHP196554	BC1 Highback Day Chair - Fabric - Tweed Charcoal	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Tweed Charcoal
CHP196556	BC1 Highback Day Chair - Fabric - Tweed Oatmeal	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Tweed Oatmeal
CHP196558	BC1 Highback Day Chair - Fabric - Tweed Denim	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Tweed Denim
CHP196560	BC1 Highback Day Chair - Fabric - Tweed	899	1	490 mm	510 mm	990 - 1180	360 - 510 mm	180 kg	Tweed Stone

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Colour
	Stone					mm			
CHP196562	BC1 Highback Day Chair - Vinyl - Ink	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Vinyl Ink
CHP196564	BC1 Highback Day Chair - Vinyl - Latte	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Vinyl Latte
CHP196568	BC1 Highback Day Chair - Vinyl - Kiwi	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Vinyl_Kiwi
CHP196570	BC1 Highback Day Chair - Vinyl - Hollywood	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Vinyl Hollywood
CHP196572	BC1 Highback Day Chair - Vinyl - Red	899	1	490 mm	510 mm	990 - 1180 mm	360 - 510 mm	180 kg	Vinyl Red
CHP196519	BC1 Highback Day Chair - King - Ink Vinyl	899	1	610 mm	610 mm	990 - 1140 mm	630 mm	220 kg	Ink Vinyl
CHP196520	BC1 Highback Day Chair - King - Latte Vinyl	899	1	610 mm	610 mm	990 - 1140 mm	630 mm	220 kg	Latte Vinyl
CHP196521	BC1 Highback Day Chair - King - Autumn Vinyl	899	1	610 mm	610 mm	990 - 1400 mm	630 mm	220 kg	Autumn Vinyl
CHP196525	BC1 Highback Day Chair - Super King - Ink Vinyl	899	1	490 mm	740 mm	990 - 1140 mm	630 mm	300 kg	Ink Vinyl
CHP196526	BC1 Highback Day Chair - Super King - Latte Vinyl	899	1	490 mm	740 mm	990 - 1140 mm	630 mm	300 kg	Latte Vinyl
CHP196527	BC1 Highback Day Chair - Super King -	879	1	490 mm	740 mm	990 - 1140	630 mm	300 kg	Autumn Vinyl

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Colour
	Autumn Vinyl					mm			
CHP196528	BC1 Highback Day Chair - Super King - Kiwi Vinyl	879	1	490 mm	740 mm	990 - 1140 mm	630 mm	300 kg	Kiwi Vinyl