

Aidacare Floorline Bed End Set



- The FL250 Bed incorporates a sophisticated range of finishing options
- From our contemporary stock selection and made to order styles through to completely bespoke, any interior design desires can be met
- Attendant Control Panel - Allows the caregiver to lock out certain functions and also ensures bed operation is not just restricted to the handset

Made to order options include:

- **Tallo** - Timber or laminate surround and recessed timber veneer inner with a straight gap detail
- **Eyre** - Timber or laminate surround and recessed timber veneer inner with a curved gap detail.
- **Como** - Timber or laminatesurround and recessed timber veneer inner.
- **Sophia** - Timber surround with upholstered inner detail on headboard and recessed timber veneer detail on footboard (please contact us for custom orders)

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	RRP	Box Quantity	Colour
BEB046200	Bed End Set - Tahoe - Sublime Teak - Single (includes: Ctrl Panel + Hardware)	609	0	-
BEB048700	Bed End Set – Como Radius – Sublime Teak – Single	495	1	Sublime Teak
BEB046210	Bed End Set - Vail - Milano Walnut - Single (includes: Ctrl Panel + Hardware)	609	0	-
BEB048701	Bed End Set – Como Radius – Milano Walnut – Single	495	1	Milano Walnut
BEB046220	Bed End Set - Tahoe - Sublime Teak - King Single (includes: Ctrl Panel +	679	0	-

Code	Display Name	RRP	Box Quantity	Colour
Hardware)				
BEB048710	Bed End Set – Como Radius – Sublime Teak – King Single	495	1	Sublime Teak
BEB046230	Bed End Set - Vail - Milano Walnut - King Single (includes: Ctrl Panel + Hardware	679	0	-
BEB048711	Bed End Set – Como Radius – Milano Walnut – King Single	495	1	Milano Walnut