

Aspire ComfiMotion Plush Pillow





- Plush Fibrefill sleep surface and Memory Foam core provides ergonomic support and comfort
- Luxurious Fibrefill Outer: Clouds of soft fibrefill provide a similar sleep experience to a traditional plush pillow
- Ventilated Memory Foam Core: Provides side and back sleepers with enhanced postural support, comfort and breathability
- Breathable Easy-Zip Cover: Encased in a removable cover for fuss-free care
- Familiar sleep surface
- Improved head and neck support
- Accomodates side and back sleepers
- Luxurious Plush Layer
- Immersive Memory Foam
- Medium Density
- Removable Cover

Warranty Information

• 1 Year

Ω	
⊘ Contract ap	proved

Code	Display Name	Total Length	Total Width	Total Height
BEP050110	Aspire ComfiMotion Plush Pillow	660 mm	480 mm	110 mm