

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
BEA007480	Alarm - Cura - Bed Monitor Plug Pack 12V DC (for 3500 Active Care Console) (replacement)	\$ 29.10	\$ 8.73	1
BEA007516	Nurse Call Test Box	\$ 120.00	\$ 36.00	1
BEA007520	Alarm - Monitor - Mono - To Suit Nurse Call Systems	\$ 184.00	\$ 55.20	2
BTC066300	Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes	\$ 210.00	\$ 63.00	1
BTC0671085	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC0671095	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC0671155	Shower Commode Elevating Legrest - Aspire - RIGHT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC069530	Shower Commode Seat - 4HC - Closed Front w/ Pan Carrier Runners - 550mm Wide	\$ 469.00	\$ 140.70	1
BTC080360	Foot Board - K Care - Shower Commode - Grey Plastic	\$ 124.00	\$ 37.20	1
BTC087636	Anti Tipppers - Invacare - for Aquatec Ocean & VIP Shower Commodes	\$ 123.00	\$ 36.90	4
BTT149700	Toilet Support Rail Toilet Roll Holder - for Devon Rail	\$ 74.00	\$ 22.20	4
BTT156500	Over Toilet Aid Toilet Seat Lid - for U Frame OTA	\$ 23.70	\$ 7.11	2
CHP213300	Hilite Chair Rear Adj Legs w/ Wheels - Vanilla/Champ	\$ 201.00	\$ 60.30	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	1
LIC356010	CP200 (Prism) Bottom Cover	\$ 43.70	\$ 13.11	3
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	1
LIC359160	Adapter for Multirall or Likorall	\$ 43.70	\$ 13.11	9
LIC360670	HandControl HB80 - Liko	\$ 374.00	\$ 112.20	1
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	1
LSS372700	Lifter Sling - Access Toileting - Kerry - Large	\$ 288.00	\$ 144.00	1
LSS373310	Lifter Sling - Invacare Hammock - Large	\$ 224.00	\$ 112.00	1
LSS374100	Sling General Purp Fab XXL (S.M.S.)	\$ 285.00	\$ 142.50	1
LSS374916	Sling - Hygiene with Head Support Polyester - Small	\$ 327.00	\$ 163.50	2
LSS376100	Sling - Supreme Hygiene - Fabric Small	\$ 327.00	\$ 163.50	2
LSS377754	Lifter Sling - Allegro - Care Head Support Mesh - Small	\$ 402.00	\$ 201.00	1
LSS379000	Lifter Sling - Seating Sling - Kerry - Large	\$ 342.00	\$ 171.00	1
LSS379900	Lifter Sling - Pivot Frame Polyester w/ Keyhole - Kerry - Small	\$ 364.00	\$ 182.00	1
LSS380900	Lifter Sling - Gen Purp Polyester - Kerry - Small	\$ 689.00	\$ 344.50	1
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
LSS394500	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Lge - 200kg	\$ 519.00	\$ 259.50	1
LSS394750	Sling - Solo Highback Sling - Liko - Med - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	1
MWC449214	Wheelchair - Start M25 - Standard - 405mm Wide (as per script)	\$ 949.00	\$ 474.50	2
MWS449205	Wheelchair Folding Light - Action 3NG - Self Propelled - 16in (405mm) - 24in solid tyres	\$ 1,329.00	\$ 664.50	1
PTA493750	HipSaver Slim Fit Veterans Kit Female X-Large	\$ 529.00	\$ 370.30	1
PTA494005	HipSaver Slim Fit Veterans Kit Male Small	\$ 529.00	\$ 370.30	1
PTA504300	Hip Protector Pants Female Size 10 (S.M.S.)	\$ 34.00	\$ 23.80	2
PTA504700	Hip Protector Pants Female Size 18 (S.M.S.)	\$ 34.00	\$ 23.80	2
PTA504900	Hip Protector Pants Female Size 22 (S.M.S.)	\$ 34.00	\$ 23.80	2
PTA505400	Hip Protector Pants Male Small (S.M.S.)	\$ 37.00	\$ 25.90	3
PTA511900	HipSaver NH Veterans Kit Male Small	\$ 529.00	\$ 370.30	1
PTA512438	HipSaver - Quick Change - Veterans Kit 2 (6-Pants, 2-Pads) - Female - Med	\$ 529.00	\$ 370.30	1
PTA512800	HipSaver - QuickChange - Starter Kit (3-Pants, 1-Pads) - Female - Med	\$ 269.00	\$ 188.30	2
PTC561274	Cushion - Jay X2 - 17 x17in	\$ 759.00	\$ 379.50	1
PTC566701	Cushion - J2 Deep Contour Cushion 24x18in	\$ 1,549.00	\$ 774.50	1
PTC571902	Cushion - JAY - Lite 14x14in	\$ 509.00	\$ 254.50	1
PTC571970	Cushion - JAY - Lite 15x16in	\$ 509.00	\$ 254.50	1
PTC594920	Cushion - ROHO Hybrid Elite 1R 16x18in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594935	Cushion - ROHO Hybrid Elite 1R 18x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC596500	Cushion Cover - ROHO Standard HP 18.5x22in	\$ 153.00	\$ 45.90	1
PTC600600	Cushion Synergy Structure 18 x 18"	\$ 772.00	\$ 386.00	1
PTM636345	ROHO Mattress Enclosure Cover	\$ 434.00	\$ 303.80	3
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	1
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTS640880	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Blk	\$ 102.00	\$ 30.60	2
PTS641310	Stockings - JOBST Relief - Thigh w/ Grip Top - 20-30mmHg - Unisex - O/Toe - Med	\$ 130.00	\$ 39.00	1
PTS641974	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Blk	\$ 107.00	\$ 32.10	1
TSE674240	Defibrillator Table - Aspire Emergency Stretchers	\$ 323.00	\$ 161.50	2
WAT907201	Avalo AX Articulating Arm with Pole Mount (needs 1x A&E Mount 12149K)	\$ 447.00	\$ 134.10	4