

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
BEA007481	Alarm - Cura - Bed Monitor Plug Pack 9V DC (replacement)	\$ 46.20	\$ 13.86	16
BEA013800	Bed Pole - Dual Prong for Single Bed	\$ 52.00	\$ 15.60	1
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2)	\$ 809.00	\$ 404.50	2
BEB041300	Bed - Alrick - P5500W - King Single 4 Section	\$ 4,219.00	\$ 1,265.70	1
BEB045900	Bed - Aidacare AC3 Single - 3 Function - CL - Beech Head & Foot Boards - 220kg	\$ 2,829.00	\$ 1,414.50	3
BEB046420	Bed End Set Custom (Sophia - Timber - No Control Panel Cut-out) - FL250/FLX - Single	\$ 879.00	\$ 263.70	1
BEC046077	Bedside Commode - Carequip - Maxi - 550mm Champagne - PC Steel - 175kg	\$ 1,669.00	\$ 834.50	1
BEC046079	Bedside Commode - Carequip - Maxi - 650mm Champagne - PC Steel - 210kg	\$ 2,219.00	\$ 1,109.50	1
BEM047704	Mattress Cover - Protect A Plas - Heavy Duty 1900x910x180mm (7in)	\$ 57.00	\$ 28.50	1
BTC062900	Shower Commode Seat - K Care - H/Made - Open 550mm - Blue	\$ 579.00	\$ 173.70	1
BTC0671105	Shower Commode Elevating Legrest - Aspire - LEFT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	2
BTC0671155	Shower Commode Elevating Legrest - Aspire - RIGHT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	2
BTC077200	Pelvic/Chest Strap - for K Care Shower Commodes - Not Fitted - Requires Rivetting	\$ 81.70	\$ 24.51	4
BTC078300	Anti Tip Bars - for K Care Shower Commodes - OLD - Compatible w Candy Cane Leg	\$ 165.60	\$ 49.68	1
BTS148095	Shower Trolley LINER, PILLOW and FITTINGS- Aspire	\$ 1,209.00	\$ 362.70	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
LSS371645	Lifter Sling - Oxford Comfort Sling Polyester - Large	\$ 404.00	\$ 202.00	1
LSS378800	Lifter Sling - Standing Lifter - Kerry - Large	\$ 297.00	\$ 148.50	1
LSS381000	Lifter Sling - Gen Purp Polyester - Kerry - Medium	\$ 609.00	\$ 304.50	1
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	1
LSS394160	Sling - SafetyVest Mod 93 - Liko - Med - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 809.00	\$ 242.70	1
LSS394190	Sling - ComfortVest Mod 95 - Liko - Med - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 799.00	\$ 239.70	1
LSS394778	Sling - Liko Universal Mod000 - Net Polyester, M	\$ 481.00	\$ 240.50	1
LSS395180	LiftSheet XL, Tall, polyester, length 185 cm (73 in)	\$ 839.00	\$ 251.70	1
MWS449600	WChair - Tilt in Space - Aspire Rehab RS - 410mm - 24in Wheels - 180kg - Gunmetal Grey	\$ 2,479.00	\$ 1,239.50	1
PTA512439	HipSaver - Quick Change - Veterans Kit 2 (6-Pants, 2-Pads) - Female - Lge	\$ 529.00	\$ 370.30	1
PTA513901	HipSaver - Nursing Home - Veterans High Compliance Pack - Female - Sml	\$ 839.00	\$ 587.30	1
PTA521697	Boot - Sheepskin Hi Temp - Pressure Care - O/Toe - XX Sml	\$ 191.00	\$ 133.70	1
PTA521898	Boot - Sheepskin Hi Temp - Short Slipper - C/Toe - X Sml	\$ 215.00	\$ 150.50	2
PTA533304M	Limb Protection - DermaSaver - Elbow Tube (2 Layers) - Med	\$ 48.70	\$ 34.09	1
PTA533313M	Limb Protection - DermaSaver - Stay-Put Heel Protector - Med	\$ 115.00	\$ 80.50	1
PTA534101	Limb/Cast Protection - Limbo - Adult Hand - Waterproof	\$ 51.00	\$ 35.70	1
PTA534102L	Limb/Cast Protection - Limbo - Adult Below Elbow - Waterproof - Lge	\$ 61.00	\$ 42.70	1
PTC549500	Cushion - Cloud - Otto Bock - 17x17in	\$ 1,279.00	\$ 639.50	1
PTC554600	Cushion - ROHO Harmony 16x18in	\$ 409.00	\$ 204.50	1
PTC573421	Cushion - JAY - Union 19x17in	\$ 649.00	\$ 324.50	1
PTC582300	JAY Fluid Pad 16x18in	\$ 929.00	\$ 464.50	1
PTC585157	Cushion - JAY - Lite 18x16	\$ 482.00	\$ 241.00	1
PTC585196	Accessories Box - JAY Cushions	\$ 1,009.00	\$ 302.70	2
PTC587400	Cushion Cover - ROHO Nexus Spirit 18x18in	\$ 179.00	\$ 89.50	1
PTC608045	Cushion - Promedicare Domino - 16 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTM636345	ROHO Mattress Enclosure Cover	\$ 434.00	\$ 303.80	1
PTS640022	JOBSTTravel Socks Black Size 2	\$ 58.00	\$ 17.40	1
PTS640024	JOBSTTravel Socks Black Size 4	\$ 58.00	\$ 17.40	2
PTS640630	Stockings - JOBST UltraSheer - Knee - 20-30mmHg - C/Toe - Lge - Blk	\$ 103.00	\$ 30.90	1
PTS640984	Stockings - JOBST UltraSheer - Waist - 30-40mmHg - P/hose - C/Toe - Med - Blk	\$ 228.00	\$ 68.40	1
PTS641210	Stockings - JOBST Relief - Knee - 30-40mmHg - Unisex - C/Toe - Lge	\$ 102.00	\$ 30.60	3
PTS641265	Stockings - Jobst Relief - Thigh w/Grip Top - 20-00mmHg - Unisex - C/Toe - Sml	\$ 130.00	\$ 39.00	1
PTS641860MRS	Stockings - GloriaMed - Unisex - Knee - 20-30mmHg - O/Toe - Med/ reg/short-Bge	\$ 101.00	\$ 30.30	2
PTS641861BE2MWL	Stockings - GloriaMed - Unisex - Knee - 20-30mmHg - C/Toe - Med/ wide/long-Bge	\$ 101.00	\$ 30.30	2
PTS641865LWL	Stockings - GloriaMed - Unisex - Knee - 30-40mmHg - O/Toe - Lge/ wide/long-Bge	\$ 103.00	\$ 30.90	1
PTS641865XLRS	Stockings - GloriaMed - Unisex - Knee - 30-40mmHg - O/Toe - XL/ reg/short-Bge	\$ 103.00	\$ 30.90	1
PTS641869RLWL	Stockings - GloriaMed - Unisex Thigh w/ Belt Top - 30-40mmHg - O/Toe - Lge /wide/long - Right - Bg	\$ 115.00	\$ 34.50	1
PTS642067BEV	Stockings - Mediven Plus - Thigh w/ Grip Top (Short) - 40-50mmHg - O/Toe - (Sz V) - Bge	\$ 230.00	\$ 69.00	1
PTS642330	Stocking - Sigvaris - Cotton Panty (ID CLASS, SIZE & COLOUR)	\$ 189.00	\$ 56.70	2
PTS647785	Stockings - Sigvaris Cotton - Knee - 23-32mmHg S Norm Nature Closed Toe	\$ 161.00	\$ 48.30	1
WAT807800	Dressing Trolley - Kerry - 500x500 - Flat Top - Stainless Steel	\$ 985.60	\$ 492.80	1
WAT907335	Drawer Divider - Tall - Grey (6 & 9in Drawers) - TRIBUTE	\$ 26.60	\$ 7.98	1