

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
BEA007214	Alarm - Invisabeam - Invisalert - Pager	\$ 839.00	\$ 251.70	3
BEA007215	Alarm - Invisabeam - Remote with Power Supply	\$ 539.00	\$ 161.70	2
BEA009320	Bed Pole - Folding - Top Mounted - 32 x 50mm	\$ 249.00	\$ 74.70	1
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	2
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	1
BEB045860	Bed - Aidacare AC3 V2 - King Single - 3 Function - Teak Bed Ends	\$ 2,950.00	\$ 1,475.00	1
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067110S	Shower Commode Elevating Legrest - Aspire - LEFT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	3
BTC087630	Chest belt	\$ 109.00	\$ 32.70	1
CHP198960	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP212700	Hilite Chair Tray Adjustable - Grey/Slate	\$ 509.00	\$ 152.70	1
CHT100400	Treatment Chair - Aspire DIALYSIS Chair - Beige Vinyl - 150kg	\$ 1,869.00	\$ 560.70	1
LIC354115	Soft Strap - P200 Weigh Scale	\$ 58.00	\$ 17.40	2
LIC355820	Limit Switch	\$ 31.30	\$ 9.39	2
LIC355840	Limit Switch Shaft	\$ 5.98	\$ 1.79	4
LIC355850	Limit Switch Roller	\$ 7.45	\$ 2.24	4
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	1
LIC359260	Extension Arm - Liko Multirall - Fits Carriage S65, Prod. No. 3136011.	\$ 194.00	\$ 58.20	1
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	1
LIC359640	Carriage with Brake for Likorall R2R for 90deg turning of motor	\$ 270.00	\$ 81.00	1
LIC360670	HandControl HB80 - Liko	\$ 374.00	\$ 112.20	1
LIC360752	FreeSpan Rail 2m - with Likorall Trolley - 250kg	\$ 1,959.00	\$ 979.50	1
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	3
PTA521913	Shoe - Sheepskin Hi Temp - DiabPro - Closed-Toe - Medium	\$ 210.00	\$ 147.00	1
PTA534102L	Limb/Cast Protection - Limbo - Adult Below Elbow - Waterproof - Lge	\$ 61.00	\$ 42.70	1
PTC573427	Cushion - JAY - Union 22 x 24in	\$ 779.00	\$ 389.50	1
PTC586900	Cushion - ROHO Nexus Spirit 16x18in	\$ 726.00	\$ 363.00	1
PTC590300	Cushion - ROHO Quadtro Select LP 11x11in (6x6cell)	\$ 1,119.00	\$ 559.50	1
PTC594921	Cushion - ROHO Hybrid Elite 1R 16x20in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTS640150	Stockings - JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Med - Blk 1pair	\$ 70.00	\$ 21.00	2
WAT907665	Avalo Writing Surface - Light Cream	\$ 173.00	\$ 51.90	1
WAT909230	Cart Accessory - Aidacare/Tribute - Sharp Container	\$ 184.00	\$ 55.20	1