

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
BEA007444	Alarm - Cura - Dual Adaptor - Black	\$ 57.00	\$ 17.10	4
BEA007473	Alarm - Cura - Monitor Soft Plastic Mounting Holder (OLD - use with 2613, 2610, 2608)	\$ 40.00	\$ 12.00	2
BEA007496	Alarm - Cura - 30 Channel Display Pad	\$ 410.00	\$ 123.00	1
BEA007514	Alarm - 6mm Extension Tele Plug - To Suit Nurse Call Systems	\$ 30.70	\$ 9.21	2
BEA007520	Alarm - Monitor - Mono - To Suit Nurse Call Systems	\$ 184.00	\$ 55.20	2
BEA007522	Alarm - Cura - Cable Connector - 6.35 to 3.5 Mono (T&S) - 2m	\$ 31.50	\$ 9.45	2
BEA007523	Alarm - Cura - Cable Connector - 6.35 Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	2
BEA007536	Alarm - Cura - Pad Connector - USB to RJ-12 Phone Plug (Wired Pin 2&5) - 2m	\$ 31.50	\$ 9.45	1
BEA007681	Alarm - Alerta - Mains Power Adaptor - Detect Receiver ONLY	\$ 48.60	\$ 14.58	5
BEA007682	Alarm - Alerta - Mains Power Adaptor - Detect Motion ONLY	\$ 48.60	\$ 14.58	5
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	8
BEA015000	Overbed/Chair Table - K Care - Tilt Dual HD - Wood/Brown	\$ 493.00	\$ 147.90	1
BEA024121	Bed End Set - CS7 wide - Drake 42" w/ Ctrl Panel & T-burg - Biltmore Cherry	\$ 789.00	\$ 236.70	4
BEA024124	Assist Bar - Invacare - for Invacare CS7 Wide Deck Bed	\$ 330.00	\$ 99.00	4
BEA024128	Assist Rail - Invacare - for Invacare CS7 Wide Deck Bed	\$ 402.00	\$ 120.60	3
BEA024131	Self Help Pole - Invacare - for Invacare All CS Beds	\$ 509.00	\$ 152.70	6
BEA029813	Mounting Bracket ONLY (NO FL250 HARDWARE) Bed Ends - Aidacare FL250 ACUTE	\$ 312.00	\$ 93.60	2
BEA029814	Moulded ABS Bed Ends - Aidacare Ward Bed/FL250 ACUTE	\$ 260.00	\$ 78.00	5
BEA826290	Bedside Cabinet - Hansel - Timber (choose stain colour)	\$ 879.00	\$ 439.50	1
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	33
BEB045105	Bed - Jomor - Cuddle Bed Configuration - 4 Function - 200kg	\$ 13,969.00	\$ 4,190.70	1
BEB045140	Bed - UC - Floorline Single	\$ 4,395.00	\$ 1,318.50	1
BEB045150	Bed - UC - Floorline Single w/ Knee Bend	\$ 5,144.90	\$ 1,543.47	3
BEB045567	Bed - Invacare - CS7-42" Wide Deck (bed ends sep)	\$ 5,659.00	\$ 2,829.50	1
BEB045569	Bed - Invacare - CS5 without head and foot board	\$ 3,709.00	\$ 1,854.50	1
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	6
BEB046295	Bed End Set Custom (Tallo - Timber (RAW) - With Control Panel Cut-out) - FL250/FLX - King Single	\$ 889.00	\$ 266.70	2
BEB046490	Bed End Set - Moulded ABS for Acute Care - Aidacare FL250 Bed	\$ 689.00	\$ 206.70	1
BEC046071	Bedside Commode - K Care - Economy - 440mm Champagne	\$ 324.00	\$ 162.00	1
BEC046130	Bedside Commode - K Care - Deluxe - 440mm Mulberry - Zinc	\$ 661.90	\$ 330.95	2
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	1
BEM046810	Mattress - Lifecomfort 1930x850x130mm (Cover Only)	\$ 257.00	\$ 128.50	1
BTC062500	Shower Commode Seat - K Care - H/Made - Closed 550mm	\$ 839.00	\$ 251.70	1
BTC062700	Shower Commode Seat - K Care - H/Made - Closed 650mm	\$ 1,029.00	\$ 308.70	3
BTC062900	Shower Commode Seat - K Care - H/Made - Open 550mm - Blue	\$ 579.00	\$ 173.70	1
BTC064100	Shower Commode Seat - Pressure-Care - K Care - H/Made - Open 600mm	\$ 1,059.00	\$ 317.70	1
BTC066300	Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes	\$ 210.00	\$ 63.00	1
BTC067103	IV Pole - Adjustable - Aspire Shower Commodes	\$ 276.00	\$ 82.80	1
BTC067109S	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	4
BTC070600SP1	SWEP/DES - Shwr Commode Transit Low w/Open Seat & Pan Carrier	\$ 1,199.00	\$ 359.70	5
BTC070900SP3	SWEP/DES - Shwr Commode SA Frest + (Paediatric) w/Closed Seat & Pan Carrier	\$ 2,519.00	\$ 755.70	3
BTC074000	Shower Commode Seat - K Care - H/Made - Open Left Side	\$ 899.00	\$ 269.70	1
BTC075300	Backrest - K Care - Curved Padded - STD Width - for K Care Shower Commodes - Fitted	\$ 310.80	\$ 93.24	1
BTC076600	Bowl Carrier - for K Care Shower Commodes - S/Steel	\$ 168.00	\$ 50.40	2
BTC077200	Pelvic/Chest Strap - for K Care Shower Commodes - Not Fitted - Requires Rivetting	\$ 81.70	\$ 24.51	1
BTC077900SP1	SWEP/DES - Shwr Commode Self Prop Low w/Open Seat & Pan Carrier	\$ 1,479.00	\$ 443.70	2
BTC077900SP2	SWEP/DES - Shwr Commode Self Prop Low w/Closed Seat & Pan Carrier	\$ 1,479.00	\$ 443.70	1
BTC079600	Contoured Headrest - for K Care Shower Commodes	\$ 999.00	\$ 299.70	1
BTC080360	Foot Board - K Care - Shower Commode - Grey Plastic	\$ 124.00	\$ 37.20	2
BTC085100	Bowl Carrier - K Care - 320mm - for K Care Shower Commodes - S/Steel	\$ 354.00	\$ 106.20	5
BTC087692	Comfort soft backrest Ocean Ergo XL	\$ 74.00	\$ 22.20	1
BTC089124	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Closed Front	\$ 346.00	\$ 103.80	1
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front	\$ 461.00	\$ 138.30	2
BTS116100	Shower Chair - Wall Mounted - Fold Away - K Care - AS Large Padded 125kg	\$ 2,358.57	\$ 707.57	1
BTS148095	Shower Trolley LINER, PILLOW and FITTINGS- Aspire	\$ 1,209.00	\$ 362.70	1
BTT146300	Over Toilet Aid - K Care - Wide 500mm 135Kg - Zinc	\$ 278.00	\$ 83.40	1
BTT148900	Urinal Non-Spill Spout	\$ 29.90	\$ 8.97	3
BTT153200	Over Toilet Aid - K Care - 445mm - Zinc	\$ 132.00	\$ 39.60	5
BTT153201	Over Toilet Aid w/Splash Guard - CareQuip - 440mm - SWL 150kg	\$ 149.00	\$ 44.70	17
BTT157118	Toilet Grab Rail - Throne Spacer Bracket - S/Steel	\$ 59.00	\$ 17.70	3
CHP192600	Chair Lowback - Westgate Bridge Chair - 450mm - pc steel	\$ 267.00	\$ 133.50	1
CHP192900	Chair Highback - Comfort Chair Deluxe - 580mm - pc steel	\$ 629.00	\$ 314.50	4
CHP193000	Chair Footrest - Adjustable - Comfort/Harbour/Westgate Chairs	\$ 184.00	\$ 92.00	6
CHP198810	Overnight Chair/Bed - Aspire SLEEPER Chair - Beige Vinyl - 150kg	\$ 3,109.00	\$ 1,554.50	4

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	6
CHP208500	Chair Lowback - K Care Utility - 440mm - Vanilla/Champ	\$ 355.00	\$ 106.50	6
CHP211100	Hilite Chair Mobile Wheels Vanilla	\$ 709.00	\$ 212.70	1
CHP211400	Hilite Chair Neckrest - Grey/Slate	\$ 436.00	\$ 130.80	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	3
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	\$ 331.00	\$ 99.30	3
CHT100100	Treatment Chair - Greiner Multiline Next IT - Navy Blue Vinyl - 220kg	\$ 6,500.00	\$ 1,950.00	1
CHT100200	Treatment Chair - Greiner Multiline Next DC - Light Grey Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	1
CHT100400	Treatment Chair - Aspire DIALYSIS Chair - Beige Vinyl - 150kg	\$ 1,869.00	\$ 560.70	2
LAC353109	Linen Bag - Half Size (OHS) - Yellow	\$ 18.70	\$ 5.61	8
LAC353780	Trolley - Kerry - Soiled Linen Trolley Single S/S w/ Foot Operated Lid	\$ 769.10	\$ 230.73	2
LAC353800	Trolley - Kerry - Soiled Linen Trolley Double No Bag & Epoxy Coated	\$ 669.00	\$ 200.70	1
LAC354050	Trolley - Kerry - Moist Linen Trolley Large Base Only, Spring-Loaded & Galvanized	\$ 719.00	\$ 215.70	1
LAC354060	Trolley - Kerry - Moist Linen Trolley Spring Base (option) - Large	\$ 728.40	\$ 218.52	1
LAC354080	Soiled Linen Trolley Bag Support Only Raised	\$ 168.00	\$ 50.40	2
LAC500555	Cleaning - Enviro Dusting Floor Mop, 350Mm,Oates	\$ 84.00	\$ 25.20	3
LAC500850	Linen Skip - Mild steel Epoxy powder coated + Plastic twin wheel - triple	\$ 619.00	\$ 185.70	3
LIC350055	PRISM TROLLEY RING GUIDE- 8 WHEEL TROLLEY - ABS (UK)	\$ 88.00	\$ 26.40	152
LIC350060	PRISM TROLLEY RING GUIDE - ABS PLASTIC (REHABASSIST)	\$ 57.00	\$ 17.10	438
LIC353550	Ceiling Hoist - Prism - Trolley for Portable Hoists - Fixed Eyelet	\$ 101.00	\$ 30.30	9
LIC354250	Ceiling Hoist - Prism - 6 Hook Carry Bar	\$ 280.00	\$ 84.00	1
LIC355480	Mascot Charger UK 9940- Prism - For Transactive XTRA	\$ 312.00	\$ 93.60	1
LIC355640	Bi-Colour Led	\$ 24.60	\$ 7.38	1
LIC355815	Battery Connector Assy - Prism - P200 Ceiling Hoist	\$ 45.50	\$ 13.65	1
LIC355820	Limit Switch	\$ 31.30	\$ 9.39	1
LIC355871	Back Cover Only - Prism - For P200	\$ 37.20	\$ 11.16	2
LIC355873	Shut Off Switch - Prism - For P200	\$ 76.00	\$ 22.80	3
LIC355882	P200 Spring Arm	\$ 6.24	\$ 1.87	4
LIC355965	CP200 CP136 HANDSET DOCK	\$ 66.00	\$ 19.80	1
LIC356005	CP200 EXTERNAL MEMBRANE SWITCH LABEL (Reval)	\$ 182.00	\$ 54.60	1
LIC356025	CP200 ON-OFF SWITCH LOOM	\$ 63.00	\$ 18.90	1
LIC356030	CP200 Motor Assembly	\$ 320.00	\$ 96.00	2
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	12
LIC359270	Extension Arm Adapter - Liko Multirall or Likorall - Fits adapter MR/LR	\$ 91.00	\$ 27.30	1
LIC359330	Spreader Bar - Liko SlingGuard 670 twin - Max 363kg (800lbs)	\$ 839.00	\$ 251.70	2
LIC359360	Hand control 2 button - Without display - LikoGuard	\$ 709.00	\$ 212.70	4
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	2
LIC359380	Standard Carriage - LikoGuard	\$ 219.00	\$ 65.70	12
LIC359390	Wall Charger - LikoGuard, AU/NZ	\$ 391.00	\$ 117.30	5
LIC359410	In Rail Charge Carriage - LikoGuard - (with Step up card and cable) - Max 363kg (800lbs)	\$ 559.00	\$ 167.70	2
LIC359420	In-Rail-Charge Charger - LikoGuard, AU	\$ 462.00	\$ 138.60	1
LIC359560	Hand Control - Likorall ES-4MT (for Lifting and Transfer Motor)	\$ 483.00	\$ 144.90	2
LIC359570	Hand Control - Likorall ES-6MTS (for Lifting, Transfer Motor and Switches)	\$ 579.00	\$ 173.70	2
LIC359800	Battery Charger - Likorall 200 (24V/0,22 A)	\$ 245.00	\$ 73.50	2
LIC360100	Rail Cover, 5 pcs a 2.0m (79 in) - Liko	\$ 89.00	\$ 26.70	3
LIC360390	Upright Support, 3,0 m (118 in), white - Liko	\$ 231.00	\$ 69.30	8
LIC360420	Turnable Bracket for Upright Support (foot, distance block, coverhats) - Liko	\$ 93.00	\$ 27.90	8
LIC360630	Side Rail Switch, electrical - White - Liko	\$ 1,749.00	\$ 524.70	2
LIC360704	Traverse Rail Carrier, Wide/Raised 100mm - Liko	\$ 809.00	\$ 242.70	1
LIC360900	UltraTwistTM Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	2
LSS363150	ArjoHuntleigh 4 PT SLING, PADDED,LGE	\$ 899.00	\$ 269.70	1
LSS367750	Sling Liko Original H/B Thin Net Poly Reinf Leg Support w Belt	\$ 442.00	\$ 221.00	1
LSS371615	Lifter Sling - Oxford - Quickfit Poly Padded - Large	\$ 459.00	\$ 229.50	1
LSS371701	Lifter Sling - Oxford - Quickfit Deluxe Padded Polyester - Lge	\$ 539.00	\$ 269.50	2
LSS374305	Sling - General Purpose - Fabric - XL	\$ 327.00	\$ 163.50	1
LSS374925	Sling - General Purpose Paediatric - Medium	\$ 261.00	\$ 130.50	1
LSS375390	Sling - Pivot Hygiene - Mesh Medium	\$ 327.00	\$ 163.50	3
LSS376300	Sling - Supreme Hygiene - Fabric Large	\$ 327.00	\$ 163.50	4
LSS378900	Lifter Sling - Standing Lifter - Kerry - w/ Buckle Large	\$ 297.00	\$ 148.50	1
LSS379000	Lifter Sling - Seating Sling - Kerry - Large	\$ 342.00	\$ 171.00	1
LSS380700	Lifter Sling - Gen Purp Mesh - Kerry - Large LR406NE3	\$ 297.00	\$ 148.50	2
LSS381000	Lifter Sling - Gen Purp Polyester - Kerry - Medium	\$ 609.00	\$ 304.50	1
LSS382400	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Medium LR408N2	\$ 376.00	\$ 188.00	3
LSS382500	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Large LR408N3	\$ 376.00	\$ 188.00	1

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
LSS383200	Lifter Sling - Hammock Mesh - Kerry - Small LR409NE1	\$ 499.00	\$ 249.50	1
LSS388800	Limb Lifter - Single Patient Use - 850mm x 140mm - 120kg per Limb (Box of 10)	\$ 209.00	\$ 62.70	10
LSS388805	Limb Lifter - Single Patient Use - Wide - 850mm x 250mm - 120kg per Limb (Box of 10)	\$ 247.00	\$ 74.10	10
LSS388810	Pannus Support - SallySling - Single Patient Use - 120kg (Box of 5)	\$ 191.00	\$ 57.30	2
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	2
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
LSS393480	Spreader Bar - Liko Universal SlingBar 350 - Max 300 kg (660 lbs).	\$ 353.00	\$ 105.90	1
LSS393540	Spreader Bar - Liko Universal SlingBar 600 - Max 300 kg (660 lbs).	\$ 490.00	\$ 147.00	1
LSS393700	Liko Adapter 12 mm - Fits Uno with serial no. from 31825, LikoLight & Viking - Max 300kg	\$ 44.90	\$ 13.47	2
LSS393800	Adapter 22 for Quick-release Hook - Liko - Fits Liko MasterLift from -89 and Likorall.	\$ 53.00	\$ 15.90	1
LSS394400	Sling - Liko Ultrasling Mod 6- Poly w/ reinforced leg suppts - short - 500kg	\$ 926.00	\$ 277.80	2
LSS394450	Sling - Liko SilhouetteSling Mod 22 - Mesh - XS - 200kg	\$ 435.00	\$ 217.50	3
LSS394460	Sling - Liko SilhouetteSling Mod 22 - Mesh - Sml - 200kg	\$ 453.00	\$ 226.50	4
LSS394500	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Lge - 200kg	\$ 519.00	\$ 259.50	1
LSS394703	Sling - Liko Original H/B MOD 200 - Plastic Coated Net - w/belt - Sml - Child	\$ 425.00	\$ 212.50	1
LSS394730	Sling - Liko ComfortSling - Mesh - XL 300kg	\$ 599.00	\$ 299.50	2
LSS394750	Sling - Solo Highback Sling - Liko - Med - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	1
MWS449205	Wheelchair Folding Light - Action 3NG - Self Propelled - 16in (405mm) - 24in solid tyres	\$ 1,329.00	\$ 664.50	1
PTA492600	HipSaver Open Bottom Veterans Kit Medium	\$ 539.00	\$ 377.30	4
PTA493400	HipSaver - SlimFit - High Compliance - Female - XXL	\$ 141.00	\$ 98.70	1
PTA493705	HipSaver - SlimFit - High Compliance - Male - X Sml	\$ 141.00	\$ 98.70	1
PTA494061	HipSaver - QuickChange - Pant Only - Male - XXL	\$ 65.00	\$ 45.50	1
PTA494065	HipSaver - Slim Fit - Pant Only - Male - XL	\$ 65.00	\$ 45.50	1
PTA505100	Hip Protector Pants Male Large (S.M.S.)	\$ 37.00	\$ 25.90	1
PTA505300	Hip Protector Pads (S.M.S.)	\$ 46.00	\$ 32.20	1
PTA509355	Heel Protector - Maxxcare Heel Pro Evolution - Medium	\$ 199.00	\$ 139.30	12
PTA509360	Heel Protector - Maxxcare Heel Pro Evolution - X Large	\$ 199.00	\$ 139.30	8
PTA511000	HipSaver - Nursing Home - Starter Kit (3-Pants, 1-Pads) - Male - Sml	\$ 269.00	\$ 188.30	1
PTA512060	HipSaver - Nursing Home - Pant Only - Male - XXL	\$ 65.00	\$ 45.50	1
PTA512520	HipSaver - Open Bottom - Pant Only - Lge	\$ 65.00	\$ 45.50	7
PTA512530	HipSaver - Open Bottom - Pant Only - Med	\$ 65.00	\$ 45.50	2
PTA512540	HipSaver - Open Bottom - Pant Only - Sml	\$ 65.00	\$ 45.50	1
PTA512947	HipSaver QC Veterans Kit Female Large	\$ 529.00	\$ 370.30	1
PTA513797	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - X Sml	\$ 170.00	\$ 119.00	1
PTA513798	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Sml	\$ 170.00	\$ 119.00	1
PTA513800	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Lge	\$ 170.00	\$ 119.00	1
PTA513801	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XL	\$ 170.00	\$ 119.00	1
PTA513802	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XXL	\$ 170.00	\$ 119.00	1
PTA513803	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - X Sml	\$ 170.00	\$ 119.00	1
PTA513804	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - Sml	\$ 170.00	\$ 119.00	1
PTA513807	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - XL	\$ 170.00	\$ 119.00	1
PTA513808	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - XXL	\$ 170.00	\$ 119.00	1
PTA520499	Boot - Sheepskin Hi Temp - Tall Slipper - Open Toe Small	\$ 276.00	\$ 193.20	1
PTA532964	DermaSaver Shin Tube Small	\$ 77.00	\$ 53.90	1
PTA533010	Limb Protection - DermaSaver - Arm Tube with Double Elbow - X Sml	\$ 93.00	\$ 65.10	1
PTA533324	Limb Protection - GelBodies - Heel & Elbow Protector - Sml	\$ 118.00	\$ 82.60	1
PTA534112M	HipSaver - Track Pant w/ Tailbone & Knee Pad - Regular - Med	\$ 208.00	\$ 145.60	1
PTC545130	Cushion - ROHO Single Compartment LP 13.25x15in (7x8cell)	\$ 1,119.00	\$ 559.50	1
PTC546400	Cushion - ROHO Single Compartment HP 16.75x9.25in (9x5cell)	\$ 1,000.00	\$ 500.00	1
PTC549900	Cushion - Cloud - Otto Bock - 19x18in	\$ 1,279.00	\$ 639.50	1
PTC551580	Cushion Vicair Adjuter 10 & incotec Cover - 16x16in	\$ 899.00	\$ 449.50	1
PTC554200	Cushion Infinity Flow Gel 18 x 16" (S.M.S.)	\$ 754.00	\$ 377.00	2
PTC556300	Cushion Cover - ROHO Waterproof LP 15x15in	\$ 179.00	\$ 53.70	1
PTC556600	Cushion Cover - ROHO Waterproof HP 16.75x18.5in	\$ 179.00	\$ 53.70	1
PTC559000	Cushion - JAY - GS 12x15in	\$ 660.00	\$ 330.00	1
PTC561271	Cushion - Jay X2 16 x18in	\$ 759.00	\$ 379.50	1
PTC561274	Cushion - Jay X2 - 17 x17in	\$ 759.00	\$ 379.50	1
PTC563100	Cushion - JAY J2 - Recline 15.5x20in	\$ 1,159.00	\$ 579.50	2
PTC569300	Cushion - JAY - Duo 16x18in	\$ 989.00	\$ 494.50	1
PTC570000	JAY Cushion Xtreme 15x16" Tritex (S.M.S.)	\$ 647.00	\$ 323.50	1
PTC571800	Cushion - JAY - Combi 16x16in	\$ 256.00	\$ 128.00	1
PTC571980	Cushion Cover - JAY - Lite 16x16in	\$ 145.00	\$ 43.50	7
PTC573800	Cushion Cover - JAY - Stretch 15.5x18in	\$ 87.00	\$ 26.10	1
PTC576200	Cushion Cover - JAY J2 - 18x20in - Recline Stretch	\$ 156.00	\$ 46.80	1
PTC576500	Cushion Cover - JAY J2 - 15.5x16in - Recline Air Exch	\$ 126.00	\$ 37.80	1
PTC577500	Cushion Cover - JAY J2 - 17x17in - Recline Stretch	\$ 126.00	\$ 37.80	1

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
PTC582100	JAY Fluid Pad Deep Contour 20x20in	\$ 1,099.00	\$ 549.50	1
PTC590400	Cushion - ROHO Quadtro Select HP 12x13in (6x7cell)	\$ 1,119.00	\$ 559.50	1
PTC594920SR	Cushion - ROHO Hybrid Elite 1R 16x18in w/HD cover Sensor Ready	\$ 1,079.00	\$ 539.50	3
PTC594921	Cushion - ROHO Hybrid Elite 1R 16x20in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594929	Cushion - ROHO Hybrid Elite 1R 17x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594935	Cushion - ROHO Hybrid Elite 1R 18x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594945	Cushion - ROHO Hybrid Elite 1R 18x20in	\$ 839.00	\$ 419.50	1
PTC594956	Cushion - ROHO Hybrid Elite 2R 15x15in	\$ 1,079.00	\$ 539.50	1
PTC600500	Cushion Synergy Structure 18 x 16"	\$ 662.00	\$ 331.00	1
PTC600600	Cushion Synergy Structure 18 x 18"	\$ 772.00	\$ 386.00	1
PTC601450	Cushion - Synergy Spectrum - Air Contour - w/ R/H valve 18 x 18in	\$ 836.20	\$ 418.10	1
PTC604100	Cushion Cover Keyhole Cushion Overslip (S.M.S.)	\$ 16.00	\$ 8.00	1
PTC609370	Cushion - Vicair - Liberty Profile Air Exchange - 18x18in	\$ 350.00	\$ 175.00	1
PTC609490	Cushion - Vicair - O2 Active 9 - 20x16in	\$ 849.00	\$ 424.50	2
PTC609800	Cushion - Vicair - Adjuster 12 O2 Incotec - 16x18in	\$ 999.00	\$ 499.50	1
PTC622510	Cushion - Pride Stealth Spectrum Cushion - Foam - 16in x 18in	\$ 467.00	\$ 233.50	1
PTM636560	Mattress Overlay - Ultimate 75 - Alt Air - 21x8cm cell - 2000x880x80mm	\$ 1,184.50	\$ 592.25	2
PTM636565	Mattress Cover - Concave Positioning Cover - Alt Air - 2000x880x100	\$ 889.00	\$ 266.70	1
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	3
PTM638015	Mattress - Aspire HYBRID - King Single - 9 Cell - 1980x1050x170mm - 180kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	2
PTS640022	JOBSTTravel Socks Black Size 2	\$ 58.00	\$ 17.40	1
PTS640024	JOBSTTravel Socks Black Size 4	\$ 58.00	\$ 17.40	1
PTS640040	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Blk	\$ 70.00	\$ 21.00	1
PTS640080	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Navy 1pair	\$ 62.00	\$ 18.60	1
PTS640270	Stockings - JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Lge - Suntan	\$ 49.00	\$ 14.70	1
PTS640315	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 15-20mmHg - C/Toe - XL - Blk	\$ 97.00	\$ 29.10	2
PTS640500	JOBST UltraShr Mater 15-20 Natural Medium	\$ 155.00	\$ 46.50	1
PTS640580	Stockings - JOBST for Men - Knee - 20-30mmHg - C/Toe - Lge - Navy	\$ 93.00	\$ 27.90	1
PTS640640	Stockings - JOBST UltraSheer - Knee - 20-30mmHg - C/Toe - XL - Blk	\$ 103.00	\$ 30.90	1
PTS640740	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 20-30mmHg - C/Toe - Sml - Natural	\$ 128.00	\$ 38.40	2
PTS640781	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Blk	\$ 207.00	\$ 62.10	1
PTS640800	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Natural	\$ 202.00	\$ 60.60	1
PTS640815	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Suntan	\$ 181.00	\$ 54.30	1
PTS640830	Stockings - JOBST for Men - Knee - 30-40mmHg - C/Toe - Sml - Blk	\$ 113.00	\$ 33.90	2
PTS640920	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Natural	\$ 102.00	\$ 30.60	2
PTS641210	Stockings - JOBST Relief - Knee - 30-40mmHg - Unisex - C/Toe - Lge	\$ 102.00	\$ 30.60	1
PTS641265	Stockings - Jobst Relief - Thigh w/Grip Top - 20-00mmHg - Unisex - C/Toe - Sml	\$ 130.00	\$ 39.00	1
PTS641335	Stockings - JOBST Relief - Waist - 20-30mmHg - P/hose - Unisex - C/Toe - XL	\$ 181.00	\$ 54.30	2
PTS641365	Stockings - JOBST Relief - Thigh w/ Belt (Chaps) Right - 20-30mmHg - Unisex - O/Toe - Lge	\$ 201.00	\$ 60.30	1
PTS641520	TED Stocking Knee Closed Toe Large Long Beige	\$ 53.00	\$ 15.90	2
PTS641550	Compress Stockings - TED - Thigh Closed Toe - Medium Long	\$ 58.00	\$ 17.40	1
PTS641632	Stockings - TED Continuing Care - Knee - Anti-embolic - C/Toe - White - Lge	\$ 63.00	\$ 18.90	2
PTS641740	Compress Stockings - TED - Knee High Small Long Beige	\$ 42.70	\$ 12.81	1
PTS641820	TED Stocking Thigh High Extra Large & Belt Long (S.M.S.)	\$ 42.00	\$ 12.60	1
PTS641860	Stockings - GloriaMed - Unisex Knee - 20-30mmHg - O/Toe - X Sml - Bge	\$ 70.00	\$ 21.00	2
PTS641957	Stockings - Venosan 4002 (Nylon) - Knee - 23-32mmHg - Unisex - O/Toe - Sml - Marokko	\$ 143.00	\$ 42.90	1
PTS641960PTMAS	Stockings - Venosan 4002 (Nylon) - Thigh / Plain - 23-32mmHg - Unisex - C/Toe - Sml - Marokko	\$ 224.00	\$ 67.20	1
PTS641974	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Blk	\$ 107.00	\$ 32.10	1
PTS641974OSBES	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Bge	\$ 107.00	\$ 32.10	1
PTS642001	Stockings - Mediven Elegance - Knee - 20-30mmHg - C/Toe - Short - (Sz I) - Bge	\$ 130.00	\$ 39.00	1
PTS642017	Stockings - Mediven Comfort - Knee (Short) - 20-30mmHg - C/Toe - (Sz I) - Bge	\$ 136.00	\$ 40.80	2
PTS648485	Stockings - Sigvaris Cotton - Thigh w/Grip Top - 18-21mmHg L Norm Nature Closed Toe	\$ 278.00	\$ 83.40	1
PTS650390	Stockings - Venosan 4001 (Nylon) - Thigh w/ Plain Top - 18-22mmHg - Unisex - O/Toe - XLge - Long - M	\$ 219.00	\$ 65.70	1
PTS650410	Stockings - Venosan 4001 (Nylon) - Waist - 18-22mmHg - Unisex - P/hose - C/Toe - XLge - Long - Mexic	\$ 255.00	\$ 76.50	1
TSE672620	Gloves Multiglide PER PAIR (S.M.S.)	\$ 37.00	\$ 11.10	2
TSE672870	Transfer - Walk Belt - Hemco - Large	\$ 45.92	\$ 13.78	17
TSE672890	Transfer - Walk Belt - Hemco - Small	\$ 45.92	\$ 13.78	17
TSE674400	Slide Sheet - SlipperySally - Single Patient Use - 2.0m x 1.45m (Box of 50)	\$ 699.00	\$ 209.70	2
TSE674410	Slide Sheet - SlipperySally - Single Patient Use (BAGGED) - 2.0m x 1.45m (Box of 50)	\$ 719.00	\$ 215.70	1
TSE674420	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1200mm wide - 454kg (Box of 5)	\$ 474.00	\$ 142.20	1
TSE674430	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1500mm wide - 454kg (Box of 5)	\$ 569.00	\$ 170.70	1
TSE674440	Air Assisted Lateral Transfer Mat Single Patient Use - 860 W x 1980 L (mm) - (Box of 10)	\$ 2,499.00	\$ 749.70	1
TSE674450	Air Assisted Lateral Transfer Mat - Air Supply 1200W	\$ 2,459.00	\$ 737.70	1
TSE674460	Air Assisted Lateral Transfer Mat - 1.7m Replacement Supply Hose	\$ 445.00	\$ 133.50	1
TSE674470	Air Assisted Lateral Transfer Mat - Transport Stand - Single Basket	\$ 1,439.00	\$ 431.70	1

AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS
BETWEEN
50-70% OFF

SALE
HURRY,
WHILST
STOCK
LASTS!

Item Number	Item Description	List Price	Now Price	Total
TSE674480	Air Assisted Lateral Transfer Mat - Disposable Waterproof Flat Sheet (Box of 50)	\$ 300.00	\$ 90.00	3
UPA200210	DDC Detergent Pump - MAXI, Optima 3	\$ 1,239.00	\$ 371.70	9
UPA200420	Temperature Monitoring Kit - Maxi & Optima 3	\$ 1,679.00	\$ 503.70	1
WAT760800	Wash Bowl - Green Autoclavable - 240mm x 110mm - 2000ml	\$ 21.40	\$ 6.42	8
WAT761400B	Bedside Cabinet - 2 Drawers + 1 Dump Drawer With Lock (Beech)	\$ 378.50	\$ 113.55	1
WAT774600	IV Stand Deluxe for wheelchairs	\$ 182.00	\$ 54.60	3
WAT805400	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Epoxy PC Frame	\$ 580.00	\$ 290.00	2
WAT805700	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Chrome Frame	\$ 1,019.00	\$ 509.50	1
WAT807800	Dressing Trolley - Kerry - 500x500 - Flat Top - Stainless Steel	\$ 985.60	\$ 492.80	2
WAT809100	Medicine Cart - Kerry - Webster Cart - 3 Drawers (48) - Beige	\$ 3,223.80	\$ 1,611.90	1
WAT816650	FlexiCart Options - Folding Shelf Heavy Duty with Flat Top	\$ 719.00	\$ 215.70	1
WAT820180	IV Pole - Aspire - 5 Star Polyurethane Base	\$ 188.00	\$ 94.00	3
WAT900735	Emergency Cart - Aidacare - 650x480x980mm - Red - 5 Drawers	\$ 3,229.00	\$ 1,614.50	3
WAT905350-O	Isolation Cart Avalo 9H - DR003 w/ No Lock & RHS Handle - YELLOW (dividers incl)	\$ 2,374.00	\$ 1,187.00	1
WAT907185	Avalo Multi-Cavity Stotage - Complete	\$ 237.00	\$ 71.10	1
WAT907250	Avalo Handle & End Panel- 10H (K) - Light Cream/Grey	\$ 425.00	\$ 127.50	3
WAT907252	Avalo Handle & End Panel- 10H (K) - Ext Red	\$ 425.00	\$ 127.50	1
WAT907254	Avalo Handle & End Panel- 10H (K) - Light Blue	\$ 425.00	\$ 127.50	1
WAT907280	Avalo End Panel only - 10H (K) - Light Cream/Grey	\$ 263.00	\$ 78.90	2
WAT907325	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Grey	\$ 2.31	\$ 0.69	4
WAT907335	Drawer Divider - Tall - Grey (6 & 9in Drawers) - TRIBUTE	\$ 26.60	\$ 7.98	7
WAT907460	Avalo Drawer Slide Assembly (K)	\$ 107.00	\$ 32.10	4
WAT907520	Avalo Keyless Pcb	\$ 390.00	\$ 117.00	1
WAT907570	Avalo Lock Tab Kit (K)	\$ 71.00	\$ 21.30	4
WAT907620	Avalo Lock Bezel Exterior Filler - Light Cream/Grey	\$ 43.90	\$ 13.17	2
WAT907665	Avalo Writing Surface - Light Cream	\$ 173.00	\$ 51.90	4
WAT907682	Avalo Bumper - Specify Colour - Ext Blue	\$ 98.00	\$ 29.40	2
WAT907716	Avalo Lock Tab 3.4 Repeat	\$ 11.90	\$ 3.57	4
WAT907717	Avalo Lock Tab 3.4 Poly	\$ 13.10	\$ 3.93	5