

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

SALE  
HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
BEA007210	Alarm - Invisabeam - Range Extender	\$ 549.00	\$ 164.70	2
BEA007411	Alarm - Cura - Plug in Transmitter- Converter to Cordless	\$ 62.60	\$ 18.78	3
BEA007414	Alarm - Cura - Premium - Bed Pad (760x255mm)	\$ 166.00	\$ 49.80	5
BEA007463	Alarm - Cura - Wireless Waterproof Call Button	\$ 327.00	\$ 98.10	1
BEA007481	Alarm - Cura - Bed Monitor Plug Pack 9V DC (replacement)	\$ 46.20	\$ 13.86	1
BEA007495	Alarm - Cura - Voice Bed Monitor	\$ 253.00	\$ 75.90	1
BEA007514	Alarm - 6mm Extension Tele Plug - To Suit Nurse Call Systems	\$ 30.70	\$ 9.21	15
BEA007522	Alarm - Cura - Cable Connector - 6.35 to 3.5 Mono (T&S) - 2m	\$ 31.50	\$ 9.45	4
BEA007523	Alarm - Cura - Cable Connector - 6.35 Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	5
BEA007528	Alarm - Cura - Cable Connector - Two Socket Input 3.5mm (Mono & Mono) to 3.5 Stereo (T&S) - 300mm	\$ 43.70	\$ 13.11	1
BEA007539	Alarm - Cura - Pad Connector - USB to 6.35mm R/ Angled Plug Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	5
BEA007681	Alarm - Alerta - Mains Power Adaptor - Detect Receiver ONLY	\$ 48.60	\$ 14.58	1
BEA007682	Alarm - Alerta - Mains Power Adaptor - Detect Motion ONLY	\$ 48.60	\$ 14.58	3
BEA007691	Alarm - Alerta - Standard Cable - Ring/Sleeve	\$ 51.00	\$ 15.30	29
BEA007692	Alarm - Alerta - Standard Cable - Tip/Ring	\$ 51.00	\$ 15.30	46
BEA007750	Bed Alarm Stand Up Nurse Call 2 Black Band Stereo Plug	\$ 809.00	\$ 242.70	7
BEA013800	Bed Pole - Dual Prong for Single Bed	\$ 52.00	\$ 15.60	3
BEA024124	Assist Bar - Invacare - for Invacare CS7 Wide Deck Bed	\$ 330.00	\$ 99.00	3
BEA024129	Assist Bar - Invacare - for Invacare CS7 Bed	\$ 323.00	\$ 96.90	5
BEA029215	Emergency Battery Back Up - Medcraft	\$ 243.00	\$ 72.90	1
BEA029814	Moulded ABS Bed Ends - Aidacare Ward Bed/FL250 ACUTE	\$ 260.00	\$ 78.00	159
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2))	\$ 809.00	\$ 404.50	2
BEA826290	Bedside Cabinet - Hansel - Timber (choose stain colour)	\$ 879.00	\$ 439.50	1
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	60
BEB045150	Bed - UC - Floorline Single w/ Knee Bend	\$ 5,144.90	\$ 1,543.47	2
BEB045152	Bed - UC - Floorline Folding - Single w/ Knee Bend	\$ 4,835.90	\$ 1,450.77	1
BEB045567	Bed - Invacare - CS7-42" Wide Deck (bed ends sep)	\$ 5,659.00	\$ 2,829.50	2
BEB045860	Bed - Aidacare AC3 V2 - King Single - 3 Function - Teak Bed Ends	\$ 2,950.00	\$ 1,475.00	1
BEC046060	Bedside Commode w/ Swing Out Bowl - K Care - 440mm Slate - Phosphate	\$ 382.90	\$ 191.45	5
BEC046120	Bedside Commode - Deluxe - 440mm Champagne - Zinc	\$ 599.00	\$ 299.50	1
BEM046736	Mattress - Aspire Lifecomfort Premier Acute - Ward Bed Single - 250kg	\$ 989.00	\$ 494.50	7
BEM047185	Mattress Overlay Action Gel Pad 690x2030mm	\$ 2,209.00	\$ 1,104.50	1
BEM047310	Convuluted Foam Mattress Overlay Cover Protect A Plus PVC 190cmx91cm	\$ 38.00	\$ 19.00	9
BEM047704	Mattress Cover - Protect A Plus - Heavy Duty 1900x910x180mm (7in)	\$ 57.00	\$ 28.50	3
BTC062400	Shower Commode Seat - K Care - H/Made - Closed 500mm	\$ 749.00	\$ 224.70	1
BTC062600	Shower Commode Seat - K Care - H/Made - Closed 600mm	\$ 999.00	\$ 299.70	1
BTC067103	IV Pole - Adjustable - Aspire Shower Commodes	\$ 276.00	\$ 82.80	2
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	2
BTC067109S	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	14
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	2
BTC071370	Shower Commode - K Care - Adjustable Height Tilt in Space - 445mm - S/Steel (Seat Sep)	\$ 4,435.24	\$ 1,330.57	1
BTC071600	Folding Safety Arms - for K Care Shower Commodes	\$ 513.90	\$ 154.17	1
BTC077200	Pelvic/Chest Strap - for K Care Shower Commodes - Not Fitted - Requires Rivetting	\$ 81.70	\$ 24.51	8
BTC078300	Anti Tip Bars - for K Care Shower Commodes - OLD - Compatible w Candy Cane Leg	\$ 165.60	\$ 49.68	1
BTC081604	Commode Seat - K Care - for K Care Bedside Commode - Padded - Slate	\$ 139.00	\$ 41.70	10
BTC083716	Shower - Self Propelling Kit for Aquatec Commode - 24in Wheel/Brakes	\$ 482.00	\$ 144.60	1
BTC087570	Commode Seat - Novis - for ETAC Mobile Shower Commode - Padded	\$ 180.00	\$ 54.00	1
BTC087608	XL Safety bar (For XL versions of the Ocean Ergo)	\$ 199.00	\$ 59.70	2
BTC089124	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Closed Front	\$ 346.00	\$ 103.80	1
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front	\$ 461.00	\$ 138.30	14
BTC089220	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 600mm - Closed Front	\$ 519.00	\$ 155.70	8
BTC089240	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 650mm - Closed Front	\$ 519.00	\$ 155.70	3
BTS118750	Mobile Shower Trolley Hose & Fitting	\$ 201.00	\$ 60.30	1
BTT157108	Toilet Grab Rail - Throne - Child Standard Rail	\$ 579.00	\$ 173.70	1
CHP190900	Chair Highback - K Care Hilite - 530mm - Grey/Slate	\$ 739.00	\$ 369.50	1
CHP191510	Chair Lowback - K Care Utility - 650mm - Grey/Slate	\$ 2,189.00	\$ 1,094.50	1
CHP192900	Chair Highback - Comfort Chair Deluxe - 580mm - pc steel	\$ 629.00	\$ 314.50	1
CHP197990	Chair Pressure Care - Seating Matters Monaco	\$ 5,305.42	\$ 2,652.71	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	1
CHP201700	Chair Footrest - Padded Large - K Care - Vanilla/Champagne	\$ 294.00	\$ 88.20	1
CHP202400	Chair Footrest - Padded Small - K Care - Vanilla/Champagne	\$ 247.00	\$ 74.10	29
CHP204500	Chair Highback - K Care Hilite - 490mm - Grey/Slate	\$ 560.00	\$ 168.00	4
CHP211400	Hilite Chair Neckrest - Grey/Slate	\$ 436.00	\$ 130.80	2
CHP213000	Hilite Chair Tray Support Tube Vanilla	\$ 265.00	\$ 79.50	3
CHP213200	Hilite Chair Rear Adj Legs w/ Wheels - Grey/Slate	\$ 201.00	\$ 60.30	2
CHP213300	Hilite Chair Rear Adj Legs w/ Wheels - Vanilla/Champ	\$ 201.00	\$ 60.30	4

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

SALE  
HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	1
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	\$ 331.00	\$ 99.30	10
CHT100210	Treatment Chair - Greiner Multiline Next DC - Navy Blue Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	2
CHT100300	Treatment Chair - Greiner Multiline XXL - Black Vinyl - 300kg	\$ 9,995.00	\$ 2,998.50	1
CHT100400	Treatment Chair - Aspire DIALYSIS Chair - Beige Vinyl - 150kg	\$ 1,869.00	\$ 560.70	16
CHT100405	Treatment Chair - Aspire DIALYSIS Chair - Navy Blue Vinyl - 150kg	\$ 1,869.00	\$ 560.70	22
FUT100190	Coffee Table - Lily Round 500 - Timber - (choose stain colour)	\$ 147.00	\$ 44.10	32
FUT100210	Coffee Table - Lily Square 500 - Timber - (choose stain colour)	\$ 147.00	\$ 44.10	17
FUT500220	Dining Table - Lily Rectangle - Timber Frame & Top - (choose stain colour)	\$ 456.00	\$ 136.80	2
FUT500250	Dining Table - Tetis Square - Timber Frame & Veneer Top - (choose stain colour)	\$ 720.00	\$ 216.00	18
FUT500260	Dining Table - Tetis Rectangle - Timber Frame & Veneer Top - (choose stain colour)	\$ 225.00	\$ 67.50	13
FUT500310	Dining Table - Peniche Rectangle - Timber Frame & Veneer Top - (choose stain colour)	\$ 559.00	\$ 167.70	3
LAC353450	Garment Rack	\$ 430.00	\$ 129.00	1
LAC353465	Coat Hangers - Plastic	\$ 1.39	\$ 0.42	49
LAC353580	Clean Linen Trolley Front Cover with Zip	\$ 433.70	\$ 130.11	10
LAC353730	Trolley - Kerry - Soiled Linen Trolley Single Epoxy Coated	\$ 224.50	\$ 67.35	5
LAC353750	Trolley - Kerry - Soiled Linen & Garbage - Single w/ Foot Operated Lid	\$ 451.70	\$ 135.51	4
LAC353770	Trolley - Kerry - Soiled Linen Trolley Single S/S	\$ 437.90	\$ 131.37	2
LAC353850	Trolley - Kerry - Soiled Linen Trolley Double S/S	\$ 739.10	\$ 221.73	1
LAC353915	Trolley - Kerry - Soiled Linen Trolley Triple S/S	\$ 1,175.30	\$ 352.59	1
LAC353980	Trolley - Kerry - Moist Linen Trolley - Small	\$ 1,969.00	\$ 590.70	2
LAC500850	Linen Skip - Mild steel Epoxy powder coated + Plastic twin wheel - triple	\$ 619.00	\$ 185.70	1
LIC350140	Ceiling Hoist - Prism Transactive TXE130D - Powered Traverse - 130kg	\$ 4,059.00	\$ 1,217.70	1
LIC350195	Ceiling Hoist - Prism Transactive TXE200DPTT (C450) - Pwr Traverse (Motor) + Pwr Turntable	\$ 4,649.00	\$ 1,394.70	3
LIC350197	Ceiling Hoist - Prism Transactive TXE200DPTXY (C450) - Pwr Traverse motor and XY	\$ 4,649.00	\$ 1,394.70	7
LIC350230	Ceiling Hoist - Prism Transactive TXE270D (C625) - Powered Traverse w/Curtain Jumper - 270kg	\$ 8,359.00	\$ 2,507.70	1
LIC350240	Ceiling Hoist - Prism Transactive TXE270D6B (C625) - Pwr Traverse w/6 Button Rtrn to Charge - 270kg	\$ 8,979.00	\$ 2,693.70	1
LIC350250	Ceiling Hoist - Prism C380M Bariatric - Manual Traverse - 380kg	\$ 6,520.00	\$ 1,956.00	1
LIC350260	Ceiling Hoist - Prism C380P Bariatric - Powered Traverse - 380kg	\$ 9,776.00	\$ 2,932.80	2
LIC351120	Ceiling Hoist - Prism - Curtain Jumper Trolley Set - Transactive Hoists	\$ 281.00	\$ 84.30	3
LIC352120	Ceiling Hoist - Prism - 2 Post Pressure Fit Gantry System	\$ 4,109.00	\$ 1,232.70	4
LIC353100	Ceiling Hoist - Prism - Track - 2.0m Rail	\$ 176.00	\$ 52.80	8
LIC353235	Ceiling Hoist - Prism - Track w/Wings - 400mm - Maintenance Opening	\$ 125.00	\$ 37.50	9
LIC353270	Ceiling Hoist - Prism - Dbl HTTrack - 5.0m Rail pearl white gloss	\$ 609.00	\$ 182.70	12
LIC353300	Ceiling Hoist - Prism - Track Curved - 15 Degree	\$ 649.00	\$ 194.70	12
LIC353320	Ceiling Hoist - Prism - Track Curved - 45 Degree	\$ 909.00	\$ 272.70	21
LIC353370	Ceiling Hoist - Prism - Charging Dock for Transactive Xtra	\$ 374.00	\$ 112.20	19
LIC353465	Ceiling Hoist - Unistrut - P1186 Angled Bracket - Zinc	\$ 12.40	\$ 3.72	28
LIC353480	Ceiling Hoist - Prism - Aluminium Box Section - Stainless Steel - 40 x 40 x 3 x 2.5m	\$ 309.00	\$ 92.70	4
LIC353535	Ceiling Hoist - Prism - Trolley for Portable Hoists - Guildman Track - No Swivel	\$ 370.00	\$ 111.00	3
LIC353550	Ceiling Hoist - Prism - Trolley for Portable Hoists - Fixed Eyelet	\$ 101.00	\$ 30.30	87
LIC353590	Ceiling Hoist - Prism - Trolley for Powered H System - XY Set	\$ 2,059.00	\$ 617.70	7
LIC353640	Ceiling Hoist - Prism - Multipoint Turntable - 16 Point - 130kg - Aluminium	\$ 2,369.00	\$ 710.70	1
LIC353650	Ceiling Hoist - Prism - Multipoint Turntable Powered - 16 Point/2 Position - 300kg - Steel	\$ 2,099.00	\$ 629.70	6
LIC354110	Weigh Scale - Marsden - for Prism P200 Ceiling Hoist	\$ 3,559.00	\$ 1,067.70	3
LIC354115	Soft Strap - P200 Weigh Scale	\$ 58.00	\$ 17.40	9
LIC354120	Ceiling Hoist - Prism - Fixed Lanyard - 30cm	\$ 104.00	\$ 31.20	6
LIC354140	Ceiling Hoist - Prism - Fixed Lanyard - 90cm	\$ 113.00	\$ 33.90	24
LIC354199	Lifting strap - Mortuary Spreader Bar - 2000x80x3 w/eyelets	\$ 99.00	\$ 29.70	5
LIC354250	Ceiling Hoist - Prism - 6 Hook Carry Bar	\$ 280.00	\$ 84.00	20
LIC354270	Ceiling Hoist - Prism - Paediatric Carry Bar	\$ 183.00	\$ 54.90	5
LIC354340	Ceiling Hoist - Prism - X-Style QRS Bariatric Carry Bar	\$ 1,339.00	\$ 401.70	1
LIC354350	Ceiling Hoist - Prism - Black Rubber Handle for QRS Reacher Bar	\$ 12.50	\$ 3.75	2
LIC355400	End Cap - 3-Hook	\$ 3.38	\$ 1.01	6
LIC355410	Easy Swivel Insert	\$ 18.00	\$ 5.40	6
LIC355420	Carry Bar Pin	\$ 0.49	\$ 0.15	2
LIC355440	Trolley Axle	\$ 12.20	\$ 3.66	2
LIC355490	Geared Wheel	\$ 23.90	\$ 7.17	4
LIC355500	Idle Wheel	\$ 8.92	\$ 2.68	4
LIC355510	Motor Mounting Plate	\$ 22.60	\$ 6.78	2
LIC355515	Angled Closing Plate With M6 Screws and Washers	\$ 22.60	\$ 6.78	5
LIC355530	Bottom Plate Assembly	\$ 539.00	\$ 161.70	1
LIC355600	2.2 X 5 Self Tapper. LCD Mount Screw	\$ 1.77	\$ 0.53	2
LIC355620	Wind Down Location Key Tx	\$ 7.67	\$ 2.30	2
LIC355640	Bi-Colour Led	\$ 24.60	\$ 7.38	1
LIC355650	Battery Series Jumper	\$ 7.03	\$ 2.11	2
LIC355670	Air Tube - 10"	\$ 8.63	\$ 2.59	5
LIC355700	Motor Mount Spacer Shim	\$ 3.40	\$ 1.02	5

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

SALE  
HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
LIC355720	6 Way Grey Rubber	\$ 8.78	\$ 2.63	4
LIC355750	Emerg. Stop Cord Assy	\$ 18.80	\$ 5.64	1
LIC355760	Slack Tape Kit	\$ 12.70	\$ 3.81	4
LIC355800	Handset Hook	\$ 14.10	\$ 4.23	8
LIC355820	Limit Switch	\$ 31.30	\$ 9.39	7
LIC355830	Carry Bar Shaft. Left Table Pin	\$ 4.01	\$ 1.20	4
LIC355860	Limit Switch Spring	\$ 19.50	\$ 5.85	4
LIC355873	Shut Off Switch - Prism - For P200	\$ 76.00	\$ 22.80	1
LIC355882	P200 Spring Arm	\$ 6.24	\$ 1.87	9
LIC355901	Brass Air Tube Inserts - Prism For P200	\$ 1.25	\$ 0.38	6
LIC355904	Hub Strap Guide - Prism For P200 - Thick White	\$ 26.50	\$ 7.95	11
LIC355905	C380 Lift Tape	\$ 799.00	\$ 239.70	2
LIC355910	CP136 Lift Tape	\$ 98.09	\$ 29.43	4
LIC355940	CP136 Australian Adapter	\$ 92.00	\$ 27.60	3
LIC355965	CP200 CP136 HANDSET DOCK	\$ 66.00	\$ 19.80	4
LIC355970	CP200 CP136 Handset Case	\$ 69.00	\$ 20.70	3
LIC355975	CP200 LIMIT SWITCH LOOM	\$ 40.00	\$ 12.00	2
LIC355980	CP200 14 TOOTH GEAR	\$ 334.00	\$ 100.20	2
LIC355985	CP200 Hub	\$ 471.00	\$ 141.30	2
LIC356000	CP200 Hoist Covers FR ABS Blue	\$ 43.70	\$ 13.11	1
LIC356005	CP200 EXTERNAL MEMBRANE SWITCH LABEL (Reval)	\$ 182.00	\$ 54.60	3
LIC356010	CP200 (Prism) Bottom Cover	\$ 43.70	\$ 13.11	1
LIC356025	CP200 ON-OFF SWITCH LOOM	\$ 63.00	\$ 18.90	3
LIC356030	CP200 Motor Assembly	\$ 320.00	\$ 96.00	1
LIC356045	Prism Track Cleaning tool - Ceiling Hoist - for PM	\$ 18.70	\$ 5.61	83
LIC357135	Fixed Lanyard 30cm	\$ 120.00	\$ 36.00	1
LIC359120	Carriage S50 - Multirall (with single hook)	\$ 252.00	\$ 75.60	11
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	1
LIC359160	Adapter for Multirall or Likorall	\$ 43.70	\$ 13.11	3
LIC359170	Extension Belt - Liko Multirall - adjustable 600-1000 mm (24-40 in)	\$ 294.00	\$ 88.20	3
LIC359200	Quick-release Hook Bolt - Multirall	\$ 219.00	\$ 65.70	2
LIC359270	Extension Arm Adapter - Liko Multirall or Likorall - Fits adapter MR/LR	\$ 91.00	\$ 27.30	3
LIC359300	Ceiling Hoist - LikoGuard L - 272kg (600lbs) SWL (no accessories)	\$ 4,919.00	\$ 1,475.70	9
LIC359320	Spreader Bar - Liko SlingGuard 450 fixed - Max 363kg (800lbs)	\$ 331.00	\$ 99.30	3
LIC359360	Hand control 2 button - Without display - LikoGuard	\$ 709.00	\$ 212.70	19
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	3
LIC359380	Standard Carriage - LikoGuard	\$ 219.00	\$ 65.70	2
LIC359410	In Rail Charge Carriage - LikoGuard - (with Step up card and cable) - Max 363kg (800lbs)	\$ 559.00	\$ 167.70	17
LIC359490	Ceiling Hoist - Likorall 200 - natural - incl hand control - 200 kg (440 lbs) SWL	\$ 2,929.00	\$ 878.70	1
LIC359535	Ceiling Hoist - Likorall 242 S - White 200Kg (440 lbs) SWL	\$ 3,119.00	\$ 935.70	2
LIC359537	Ceiling Hoist - Likorall 243 ES - White 230Kg (507 lbs) SWL	\$ 3,669.00	\$ 1,100.70	2
LIC359560	Hand Control - Likorall ES-4MT (for Lifting and Transfer Motor)	\$ 483.00	\$ 144.90	4
LIC359570	Hand Control - Likorall ES-6MTS (for Lifting, Transfer Motor and Switches)	\$ 579.00	\$ 173.70	10
LIC359610	Quick-release Carriage Likorall - Fits Quick-release Adapter, Prod. No. 3126009	\$ 159.00	\$ 47.70	5
LIC359660	Carriage Adapter - Likorall for Quick-release Hook Fits Carriage S50, Prod. No. 3136010	\$ 125.00	\$ 37.50	9
LIC359760	Extension Flex Multi 1300 - Extra long charger cable - Traverse Carrier Multi and MultiStation	\$ 130.00	\$ 39.00	13
LIC359800	Battery Charger - Likorall 200 (24V/0,22 A)	\$ 245.00	\$ 73.50	4
LIC359885	3.4m Straight Rail H70 - White - Liko	\$ 443.00	\$ 132.90	1
LIC359940	4m Embedded Straight Rail H70 - White - Liko	\$ 539.00	\$ 161.70	5
LIC359950	90deg Embedded Rail Curve H70 - White, 1m/1m	\$ 649.00	\$ 194.70	4
LIC359960	45deg Embedded Rail Curve H70 - White, 0.6m/0.6m	\$ 483.00	\$ 144.90	2
LIC360100	Rail Cover, 5 pcs a 2.0m (79 in) - Liko	\$ 89.00	\$ 26.70	28
LIC360160	Pendant adjustable 90-130 mm (4-5 in) - White - Liko	\$ 124.00	\$ 37.20	13
LIC360320	Pendant Cover dia 50, 2-pack (2.0 in.) - Liko	\$ 13.80	\$ 4.14	584
LIC360330	Pendant Washer 1 mm (0.04 in.), 50 pcs - Liko	\$ 12.50	\$ 3.75	4
LIC360350	Slimline Upright Support, L 3 m (118 in) - White - Liko	\$ 193.00	\$ 57.90	20
LIC360370	Slimline Foot, high (add on assembly parts Prod. No. 3101661) - Liko	\$ 408.00	\$ 122.40	6
LIC360390	Upright Support, 3,0 m (118 in), white - Liko	\$ 231.00	\$ 69.30	13
LIC360410	Upright Support Distance Aluspacer - Liko	\$ 231.00	\$ 69.30	2
LIC360420	Turnable Bracket for Upright Support (foot, distance block, coverhats) - Liko	\$ 93.00	\$ 27.90	6
LIC360440	Upright Support Steel, 2.6 m (102 in), white - Liko	\$ 322.00	\$ 96.60	21
LIC360480	Wall Bracket, parallel - Liko	\$ 115.00	\$ 34.50	4
LIC360510	Assembly Parts - Suspended Rail System - Liko - Bayonet Rail 1.0m, white Fits rails H140/160/180	\$ 155.00	\$ 46.50	6
LIC360520	Assembly Set reinforcement rail 50 - Liko	\$ 115.00	\$ 34.50	8
LIC360540	Reinforcement Rail 30, 1,45 m, white (58 in.) - Liko	\$ 356.00	\$ 106.80	6
LIC360550	Reinforcement Rail 50, 1,95 m, white (77 in.) - Liko	\$ 959.00	\$ 287.70	2
LIC360630	Side Rail Switch, electrical - White - Liko	\$ 1,749.00	\$ 524.70	2
LIC360640	Turntable, mechanical - White - Liko	\$ 2,309.00	\$ 692.70	2

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
LIC360650	Turntable, electrical - White - Liko	\$ 3,309.00	\$ 992.70	3
LIC360680	Extension cable HB80 - Liko	\$ 287.00	\$ 86.10	1
LIC360707	Traverse Rail Carrier, extra wide - Liko	\$ 1,619.00	\$ 485.70	2
LIC360710	Traverse Carrier Multi, standard - Incl. 1 pce Traverse Rail Carrier, standard. Fits Likorall ES.	\$ 1,329.00	\$ 398.70	4
LIC360720	Traverse Carrier Multi, wide - Incl. 1 pce Traverse Rail Carrier, wide. Fits Likorall ES.	\$ 1,329.00	\$ 398.70	3
LIC360752	FreeSpan Rail 2m - with Likorall Trolley - 250kg	\$ 1,959.00	\$ 979.50	2
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	1
LIC360770	FreeSpan Rail 3m - with Multirall Trolley - 250kg	\$ 2,279.00	\$ 1,139.50	1
LIC360870	Ultra Control Unit - incl, hand control - Fits Likorall ES and rails	\$ 1,629.00	\$ 488.70	1
LIC360900	UltraTwist™ Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	3
LIC360960	Travers Rail Carrier, E-System (Wide) - Liko	\$ 1,269.00	\$ 380.70	1
LIC360980	Curtain Rail T, 2.5m - T-System - Liko	\$ 569.00	\$ 170.70	6
LIC361050	Pendant T, adjustable 1100-1350mm - T-System - Liko	\$ 331.00	\$ 99.30	3
LIC361200	Electrical Motor for Turntable - Liko	\$ 699.00	\$ 209.70	1
LIC361210	Wedge (key for motor shaft) for Turntable - Liko	\$ 6.24	\$ 1.87	1
LIC361220	Friction Ring for Turntable - Liko	\$ 71.00	\$ 21.30	1
LSS361100	Sling Quickfit Deluxe Small (S.M.S.)	\$ 468.00	\$ 140.40	1
LSS363162	ArjoHuntleigh 4PT TOILETING SLING H/S Extra Large	\$ 869.00	\$ 260.70	1
LSS365550	Travel Bag for LikoLight Lifter	\$ 861.50	\$ 258.45	4
LSS367750	Sling Liko Original H/B Thin Net Poly Reinf Leg Support w Belt	\$ 442.00	\$ 221.00	3
LSS367864	Sling - Liko Original H/B MOD 200 - Polyester With Belt - Sml	\$ 443.00	\$ 221.50	1
LSS371629	Lifter Sling - Oxford - Quickfit Poly Padded - Medium	\$ 229.00	\$ 114.50	5
LSS371698	Lifter Sling - Oxford - Quickfit Deluxe Poly - Small	\$ 217.00	\$ 108.50	1
LSS373400	Sling General Purpose Fab Small	\$ 327.00	\$ 163.50	1
LSS373495	Sling - General Purp Mesh - Small	\$ 327.00	\$ 163.50	1
LSS374590	Sling - Pivot - Mesh - XS	\$ 311.00	\$ 155.50	1
LSS374900	Sling - Pivot - Mesh - XL	\$ 261.00	\$ 130.50	1
LSS374920	Sling - General Purpose Paediatric - Small	\$ 261.00	\$ 130.50	1
LSS376100	Sling - Supreme Hygiene - Fabric Small	\$ 327.00	\$ 163.50	1
LSS376405	Sling - Standup - Roze - Buckle - Loops - Medium	\$ 261.00	\$ 130.50	1
LSS376415	Sling - Standup - Roze - Velcro - Loops - Universal Size	\$ 261.00	\$ 130.50	1
LSS377390	Lifter Pivot Frame - 4 point - to suit Stature Hoist	\$ 559.00	\$ 167.70	1
LSS377415	Lifter Yoke to suit Alto Hoist	\$ 205.00	\$ 61.50	1
LSS377769	Lifter Sling - Allegro - Head Support Pivot - Mesh - Large	\$ 444.00	\$ 222.00	3
LSS378451	Lifter Leg Assembly to Suit Kerry KH405F	\$ 366.00	\$ 109.80	1
LSS378500	Hoist Leg Strap & Clip On Strap	\$ 156.10	\$ 46.83	2
LSS386202	Lifter - Invacare - 2 Point Spreader Bar - Jasmine	\$ 262.00	\$ 78.60	2
LSS386203	Lifter - Invacare - Adaption Kit - Jasmine	\$ 326.00	\$ 97.80	1
LSS387550	Scales Platform Electronic 360kg, hand rail on one side w/access ramps	\$ 6,249.00	\$ 1,874.70	1
LSS388800	Limb Lifter - Single Patient Use - 850mm x 140mm - 120kg per Limb (Box of 10)	\$ 209.00	\$ 62.70	1
LSS388805	Limb Lifter - Single Patient Use - Wide - 850mm x 250mm - 120kg per Limb (Box of 10)	\$ 247.00	\$ 74.10	2
LSS388810	Pannus Support - SallySling - Single Patient Use - 120kg (Box of 5)	\$ 191.00	\$ 57.30	1
LSS390916	Lifter Sling - Aspire - SML - Deluxe General Purpose w/Head Supp - CLIPS - Mesh - 200kg	\$ 423.00	\$ 211.50	8
LSS390940	Lifter Sling - Aspire - SML - Deluxe Hygiene/Access w/Head Supp - CLIPS - Poly - 200kg	\$ 480.00	\$ 240.00	18
LSS390960	Lifter Sling - Aspire - SML - Low Friction General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	20
LSS390965	Lifter Sling - Aspire - MED - All Day General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	1
LSS391500	Lifter - Patient Lifter - Essentials 180 - Electric Leg - 180kg	\$ 2,489.00	\$ 1,244.50	1
LSS393130	Low Base Kit - Golvo 7007/8008/9000	\$ 1,579.00	\$ 473.70	3
LSS393210	Leg Protector for Viking L	\$ 59.00	\$ 17.70	5
LSS393310	Sabina II Heel Support (Fits Sabina II and Sabina 200)	\$ 359.00	\$ 107.70	5
LSS393350	Sabina Calf Strap	\$ 151.00	\$ 45.30	2
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	6
LSS393410	Wall Mounted Battery Charger, AU/NZ for CBL Battery.	\$ 679.00	\$ 203.70	5
LSS393455	Wheel Set 100/75, (lower clearance) - Fits Viking L & XL	\$ 250.00	\$ 75.00	6
LSS393480	Spreader Bar - Liko Universal SlingBar 350 - Max 300 kg (660 lbs).	\$ 353.00	\$ 105.90	8
LSS393490	Spreader Bar - Liko Universal SlingBar 350 with Quick-release Hook - Max 300 kg (660 lbs).	\$ 431.00	\$ 129.30	56
LSS393540	Spreader Bar - Liko Universal SlingBar 600 - Max 300 kg (660 lbs).	\$ 490.00	\$ 147.00	4
LSS393590	Sling Cross-bar 450 with Quick-release Hook - Liko - Max 300kg	\$ 879.00	\$ 263.70	4
LSS393610	Sling Cross-bar 670 with Quick-release Hook - Liko - Max 300kg	\$ 1,019.00	\$ 305.70	2
LSS393800	Adapter 22 for Quick-release Hook - Liko - Fits Liko MasterLift from -89 and Likorall.	\$ 53.00	\$ 15.90	4
LSS394120	Sling - SupportVest w/ padding Mod 91 - Liko - Sml - w/2 Loop Clips Fits Sabina II 350 SlingBar.	\$ 709.00	\$ 212.70	2
LSS394150	Sling - SupportVest w/ padding Mod 91 - Liko - XL - w/2 Loop Clips Fits Sabina II 350 SlingBar.	\$ 749.00	\$ 224.70	2
LSS394160	Sling - SafetyVest Mod 93 - Liko - Med - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 809.00	\$ 242.70	5
LSS394190	Sling - ComfortVest Mod 95 - Liko - Med - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 799.00	\$ 239.70	10
LSS394200	Sling - ComfortVest Mod 95 - Liko - Lge - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 829.00	\$ 248.70	3
LSS394210	Sling - ComfortVest Mod 95 - Liko - XL - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 849.00	\$ 254.70	4
LSS394220	Sling - HygieneVest Highback Mod 55 - Liko - Poly w/ reinforced leg supports - Sml - 200kg	\$ 1,029.00	\$ 308.70	1
LSS394240	Sling - HygieneVest Highback Mod 55 - Liko - Poly w/ reinforced leg supports - Lge - 200kg	\$ 1,129.00	\$ 338.70	2

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
LSS394300	Sling (Walking) - MasterVest - Front Close Mod 64 - Liko - Poly - XXL - 200kg	\$ 929.00	\$ 278.70	6
LSS394395	Solo Cover For Ultra LiftPants Mod 92 , XXL/XXXL, 2x10pcs Nonwoven Material, Not Washable	\$ 669.00	\$ 200.70	21
LSS394450	Sling - Liko SilhouetteSling Mod 22 - Mesh - XS - 200kg	\$ 435.00	\$ 217.50	14
LSS394460	Sling - Liko SilhouetteSling Mod 22 - Mesh - Sml - 200kg	\$ 453.00	\$ 226.50	10
LSS394480	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Sml - 200kg	\$ 482.00	\$ 241.00	6
LSS394500	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Lge - 200kg	\$ 519.00	\$ 259.50	2
LSS394550	Sling - Liko Amputee Highback - Poly - Reinforced leg supports - XL - 300kg	\$ 579.00	\$ 289.50	1
LSS394660	Sling - Liko Original H/B MOD 200 - Net Polyester - Reinforced Leg Supports - w/belt - XS - Child	\$ 519.00	\$ 259.50	9
LSS394674	Sling - Liko Original H/B MOD 200/210 - Net Polyester - Med Slim (Slim Adult)	\$ 432.00	\$ 216.00	10
LSS394675	Sling - Liko Original H/B MOD 200/210 - Net Polyester - Reinforced Leg Supports - Med Slim (Slim Adu	\$ 559.00	\$ 279.50	13
LSS394708	Sling - Liko Original H/B - MOD200/210 - Plastic Coated Net (Bath) - XL (max 500Kg)	\$ 539.00	\$ 269.50	2
LSS394730	Sling - Liko ComfortSling - Mesh - XL 300kg	\$ 599.00	\$ 299.50	13
LSS394736	Sling - Liko - ComfortSling Plus 300 - Mesh - XL	\$ 519.00	\$ 259.50	2
LSS394740	Sling - Solo Highback Sling - Liko - Med (Slim) - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	100
LSS394750	Sling - Solo Highback Sling - Liko - Med - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	9
LSS394760	Sling - Solo Highback Sling - Liko - Lge - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	9
LSS394770	Sling - Solo Highback Sling - Liko - XLge - 2x4 pcs - No Wash - 500kg	\$ 1,379.00	\$ 689.50	56
LSS394772	Sling - Liko Universal Mod000 - net polyester, reinforced leg supports, L	\$ 579.00	\$ 289.50	2
LSS394778	Sling - Liko Universal Mod000 - Net Polyester, M	\$ 481.00	\$ 240.50	7
LSS395060	Spatula Mod 600 IC	\$ 189.00	\$ 56.70	1
LSS395110	SideRails 600 IC - Liko - 250kg	\$ 779.00	\$ 233.70	1
LSS395220	Stretch Leveller - Fits Mod 600 IC and OctoStretch - 300kg	\$ 839.00	\$ 251.70	3
LSS395370	FixStrap 135cm (10pc set) - Liko	\$ 179.00	\$ 53.70	10
MWS441400	Wheelchair Folding - Aspire ASSIST - 350 Silver	\$ 449.00	\$ 224.50	5
MWS448825	Glide Series 3 Heavy Duty Wheelchair 49cm	\$ 2,519.00	\$ 1,259.50	1
MWS449091	Wheelchair Breezy Basix2 14in - 24in Solid Wheels - RED	\$ 819.00	\$ 409.50	1
PTA493750	HipSaver Slim Fit Veterans Kit Female X-Large	\$ 529.00	\$ 370.30	1
PTA493800	HipSaver Slim Fit Veterans Kit Female Large	\$ 529.00	\$ 370.30	1
PTA494154	HipSaver - QuickChange - Pant Only - Female - Lge	\$ 65.00	\$ 45.50	1
PTA512428	HipSaver - SlimFit - Veterans Kit 2 (6-Pants, 2-Pads) - Female - XL	\$ 529.00	\$ 370.30	3
PTA512438	HipSaver - Quick Change - Veterans Kit 2 (6-Pants, 2-Pads) - Female - Med	\$ 529.00	\$ 370.30	1
PTA513797	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - X Sml	\$ 170.00	\$ 119.00	1
PTA513801	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XL	\$ 170.00	\$ 119.00	2
PTA513802	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XXL	\$ 170.00	\$ 119.00	1
PTA513803	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - X Sml	\$ 170.00	\$ 119.00	1
PTA513808	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - XXL	\$ 170.00	\$ 119.00	1
PTA513823	HipSaver - SlimFit - Veterans High Compliance Pack - Female - Lge	\$ 839.00	\$ 587.30	1
PTA513883	HipSaver - SlimFit - Veterans High Compliance Pack - Male - Lge	\$ 839.00	\$ 587.30	1
PTA532966	DermaSaver Shin Tube Large	\$ 77.00	\$ 53.90	8
PTA533318L	Limb Protection - DermaSaver - Leg Tube with Double Knee - Lge	\$ 132.00	\$ 92.40	1
PTA533324	Limb Protection - GelBodies - Heel & Elbow Protector - Sml	\$ 118.00	\$ 82.60	1
PTA534106	Limb/Cast Protection - Limbo - Adult Knee Injury - WaterproofM86	\$ 51.00	\$ 35.70	1
PTA534111L	HipSaver - Track Pant w/ Knee Pad - Regular - Lge	\$ 181.00	\$ 126.70	1
PTC545377	Cushion - ROHO Single Compartment LP 15x11.5in (8x6cell)	\$ 1,000.00	\$ 500.00	1
PTC556400	Cushion Cover - ROHO Waterproof HP 15x16.75in	\$ 179.00	\$ 53.70	3
PTC558600	Cushion - JAY - GS 14x17in w/ Incontinence Cover	\$ 1,119.00	\$ 559.50	1
PTC561000	Cushion - JAY - Combi w/ Lumbar 16x18in	\$ 301.00	\$ 150.50	2
PTC563500	Cushion - JAY J2 - 14x16in w/ Incontinence Cover	\$ 1,169.00	\$ 584.50	1
PTC567200	Cushion - JAY J2 - Plus Bariatric 22x20in	\$ 1,609.00	\$ 804.50	1
PTC567800	Cushion - JAY - Care 18x15.5in	\$ 949.00	\$ 474.50	2
PTC571901	Cushion - JAY - Lite 14x16in	\$ 509.00	\$ 254.50	1
PTC573400	Cushion - JAY - Combi Soft 12x12in	\$ 212.00	\$ 106.00	1
PTC580004	Cushion - JAY - Easy 14x16in Curved Bottom	\$ 779.00	\$ 389.50	1
PTC585157	Cushion - JAY - Lite 18x16	\$ 482.00	\$ 241.00	1
PTC590400	Cushion - ROHO Quadro Select HP 12x13in (6x7cell)	\$ 1,119.00	\$ 559.50	1
PTC591600	Cushion - ROHO Quadro Select HP 15x13in (8x7cell)	\$ 1,119.00	\$ 559.50	1
PTC594945	Cushion - ROHO Hybrid Elite 1R 18x20in	\$ 839.00	\$ 419.50	1
PTC594946	Cushion - ROHO Hybrid Elite 1R 22x18in Sensor Ready	\$ 1,239.00	\$ 619.50	1
PTC608040	Cushion - Promedicare Domino - 16 x 16in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608050	Cushion - Promedicare Domino - 18 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608056	Cushion - Promedicare Domino - 19 x 18in w/ 3D cover	\$ 1,619.00	\$ 809.50	1
PTC608322	Cushion Cover - Tanto - 3D - 16 x 16in	\$ 294.00	\$ 88.20	1
PTM636345	ROHO Mattress Enclosure Cover	\$ 434.00	\$ 303.80	1
PTM636360	Mattress Overlay - ROHO Levelling Pad - 34inx 80in	\$ 709.00	\$ 354.50	1
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	3
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTM638025	Mattress - Aspire HYBRID ACUTE - K/Single - 9 Cell - 1980x1050x170mm - 230kg (Pump Sep: PTM638030)	\$ 2,155.00	\$ 1,077.50	1
PT5640040	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Blk	\$ 70.00	\$ 21.00	1

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
PTS641860M	Stockings - GloriaMed - Unisex Knee - 20-30mmHg - O/Toe - Med - Bge	\$ 70.00	\$ 21.00	1
PTS641873MWS	Stockings - GloriaMed - Unisex (Cotton) - Knee - 20-30mmHg - O/Toe - Med/ wide/short - Bge	\$ 81.90	\$ 24.57	4
PTS641881	Stockings - GloriaMed - Lady Gloria - Thigh w/ Grip Top - 15-20mmHg - C/Toe - (Sz 1) - Bge	\$ 74.00	\$ 22.20	2
PTS641881BL4	Stockings - GloriaMed - Lady Gloria - Thigh w/ Grip Top - 15-20mmHg - C/Toe - (Sz 4) - Blk	\$ 74.00	\$ 22.20	2
PTS641978OLBEXL	Stockings - Venosan 6002 (Nyl/Lyc) - Thigh / Plain - 23-32mmHg - Unisex - O/Toe - Long - XL - Bge	\$ 188.00	\$ 56.40	1
PTS641978OLBLXL	Stockings - Venosan 6002 (Nyl/Lyc) - Thigh / Plain - 23-32mmHg - Unisex - O/Toe - Long - XL - Blk	\$ 188.00	\$ 56.40	1
PTS641997BEM	Stockings - Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Med - Bge	\$ 61.00	\$ 18.30	1
PTS642019BEIV	Stockings - Mediven Comfort - Knee (Short) - 20-30mmHg - O/Toe - (Sz IV) - Bge	\$ 136.00	\$ 40.80	1
PTS642103BEM	Stockings - DUOMed - Knee - 20-30mmHg - C/Toe - Med - Bge	\$ 78.00	\$ 23.40	1
PTS642104	Stockings - DUOMed - Knee - 30-40mmHg - O/Toe - Sml - Blk	\$ 78.00	\$ 23.40	1
PTS6426155	Solaris - ExoStrong - Thigh - 20 to 30mmHg - Short - Small	\$ 285.00	\$ 85.50	1
TSE672200	Transfer - Board - Slide Board - Light 1500Lx520W	\$ 529.00	\$ 264.50	3
TSE672526-S	Evacuation Sheet - Aidacare Lifecomfort - Single - 1900mm x 720mm - 300kg	\$ 150.00	\$ 75.00	2
TSE672527	Transfer Tube	\$ 36.10	\$ 10.83	1
TSE672770	Transfer - Patient Turntable - 390mm (Diameter) (18mm Thick)	\$ 229.00	\$ 68.70	1
TSE672870	Transfer - Walk Belt - Hemco - Large	\$ 45.92	\$ 13.78	1
TSE672890	Transfer - Walk Belt - Hemco - Small	\$ 45.92	\$ 13.78	8
TSE674220	Hospital Stretcher - Aspire Advanced Emergency Stretcher - WIDE - Hydraulic - 225kg	\$ 8,899.00	\$ 4,449.50	20
TSE674240	Defibrillator Table - Aspire Emergency Stretchers	\$ 323.00	\$ 161.50	3
TSE674400	Slide Sheet - SlipperySally - Single Patient Use - 2.0m x 1.45m (Box of 50)	\$ 699.00	\$ 209.70	2
TSE674410	Slide Sheet - SlipperySally - Single Patient Use (BAGGED) - 2.0m x 1.45m (Box of 50)	\$ 719.00	\$ 215.70	2
TSE674420	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1200mm wide - 454kg (Box of 5)	\$ 474.00	\$ 142.20	1
TSE674430	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1500mm wide - 454kg (Box of 5)	\$ 569.00	\$ 170.70	1
TSE674440	Air Assisted Lateral Transfer Mat Single Patient Use - 860 W x 1980 L (mm) - (Box of 10)	\$ 2,499.00	\$ 749.70	2
TSE674450	Air Assisted Lateral Transfer Mat - Air Supply 1200W	\$ 2,459.00	\$ 737.70	2
TSE674460	Air Assisted Lateral Transfer Mat - 1.7m Replacement Supply Hose	\$ 445.00	\$ 133.50	1
TSE674470	Air Assisted Lateral Transfer Mat - Transport Stand - Single Basket	\$ 1,439.00	\$ 431.70	2
TSE674480	Air Assisted Lateral Transfer Mat - Disposable Waterproof Flat Sheet (Box of 50)	\$ 300.00	\$ 90.00	1
UPA200100	DDC Detergent Basket - Lockable	\$ 403.00	\$ 120.90	7
UPA200210	DDC Detergent Pump - MAXI, Optima 3	\$ 1,239.00	\$ 371.70	9
UPA200230	DDC Detergent Pump - Incomatic	\$ 1,239.00	\$ 371.70	5
UPA200240	DDC - Optima 2 / Midi - Utensil tray + sides 2 pc set	\$ 142.00	\$ 42.60	7
UPA200276	Bed Pan Rack w/ Drip Tray (4 Pan / 4 Bottle) - S/Steel	\$ 1,059.00	\$ 317.70	3
UPA200285	Utility Wall Shelf - 72 x 37 x 33cm deep - S/Steel	\$ 579.00	\$ 173.70	1
UPA200300	DDC - Optima 2/Midi Chemica Monitoring System	\$ 1,819.00	\$ 545.70	2
UPA200310	DDC - Optima 3/Maxi Chemical Monitoring System	\$ 1,819.00	\$ 545.70	1
UPA200420	Temperature Monitoring Kit - Maxi & Optima 3	\$ 1,679.00	\$ 503.70	1
UPA200430	Temperature Monitoring Kit - Maxi & Optima 2	\$ 1,679.00	\$ 503.70	4
UPP100110	Plasticware - Bed Pan Support for (91-01-005)	\$ 57.00	\$ 17.10	60
UPP100230D	Pulp - Kidney Dish - 700ml	\$ 137.00	\$ 41.10	5
UPP100244	Pulp - Fluid Absorber - 6g Sachet	\$ 170.00	\$ 51.00	13
UPP100245	Pulp - Male Urine Bottle - 900ml - Square	\$ 149.00	\$ 44.70	27
WAT773600	FlexiCart Options - Folding Shelf with Flat Top	\$ 390.70	\$ 117.21	1
WAT774600	IV Stand Deluxe for wheelchairs	\$ 182.00	\$ 54.60	1
WAT807300	Dressing Trolley - Kerry - 800x500 - 1 S/S Drawer - Chrome Frame	\$ 967.00	\$ 483.50	3
WAT807400	Dressing Trolley - Kerry - 800x500 - 2 S/S Drawers - Chrome Frame	\$ 2,002.70	\$ 1,001.35	1
WAT809100	Medicine Cart - Kerry - Webster Cart - 3 Drawers (48) - Beige	\$ 3,223.80	\$ 1,611.90	1
WAT812200	Emergency Cart - Kerry - FlexiCart - 5 Drawer - Oxy Carrier / IV Pole	\$ 6,329.00	\$ 3,164.50	2
WAT814900	FlexiCart Options - Chart Holder - Standard A4	\$ 233.00	\$ 69.90	2
WAT817500	FlexiCart Options - Lockrail	\$ 386.00	\$ 115.80	1
WAT818500	FlexiCart Options - Waste Bin (No Bracket)	\$ 105.00	\$ 31.50	2
WAT819010	Catheter Floor Stand - 230 x 500mm	\$ 62.00	\$ 18.60	6
WAT820100	IV Pole PH - Kerry	\$ 155.00	\$ 46.50	2
WAT820180	IV Pole - Aspire - 5 Star Polyurethane Base	\$ 188.00	\$ 94.00	2
WAT821560	Exam Couch - Aidacare - 3 Section Electric Hi/Lo w/Short Head - Navy Blue	\$ 2,619.00	\$ 1,309.50	1
WAT822620	Examination Couch - Electric Adjustable Height - 2055x650x585-935mm - 175kg	\$ 2,489.00	\$ 1,244.50	1
WAT822650	Examination Table - Aspire Baby Change Cabinet - 1100x550x970mm - 50kg	\$ 1,195.00	\$ 597.50	1
WAT822900	Post Op Step Single - Kerry - Grey Large	\$ 120.00	\$ 60.00	1
WAT823900	Food Dispenser Trolley Tarp Cover - Kerry - 16 Tray	\$ 973.40	\$ 486.70	4
WAT900605	Instrument Trolley - Aidacare - S/Steel - 800x500x900 - 2 Shelves	\$ 869.00	\$ 434.50	1
WAT900653	Dressing Trolley - S/Steel - 800 x 500 x 900mm - 2 Drawer	\$ 1,839.00	\$ 919.50	1
WAT900660	Dressing Trolley - Aidacare - S/S - 800x500x900 - 4 Drawer	\$ 1,799.00	\$ 899.50	1
WAT900710	Treatment Cart - Aidacare - 650x480x980mm - Blue - 2 Drawers	\$ 1,989.00	\$ 994.50	1
WAT900715	Treatment Cart - Aidacare - 650x480x980mm - Blue - 3 Drawers	\$ 2,109.00	\$ 1,054.50	1
WAT900735	Emergency Cart - Aidacare - 650x480x980mm - Red - 5 Drawers	\$ 3,229.00	\$ 1,614.50	1
WAT900740	Chart Trolley - Aidacare - 650x480x900mm - Blue - No Drawers	\$ 2,359.00	\$ 707.70	2
WAT903195	Rack - 1 bay - stainless steel + rollers	\$ 2,980.00	\$ 894.00	1
WAT905180	Med Cart TRIBUTE 110 - 10H DR131 - Digital-Lock & RHS Handle - BLUE (dividers incl)	\$ 2,879.00	\$ 1,439.50	2

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

SALE  
HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	List Price	Now Price	Total
WAT905185	Med Cart TRIBUTE NS9 - NARROW - L580 x W540 x H1000mm - Key Lock - CREAM (dividers incl)	\$ 2,285.78	\$ 1,142.89	1
WAT905245	Med Cart TRIBUTE NS9 - NARROW - L580 x W540 x H1000mm - Digital-Lock - CREAM (dividers incl)	\$ 2,285.78	\$ 1,142.89	1
WAT905270	Med Cart TRIBUTE 110 - 10H DR131 - Digital-Lock & RHS Handle - CREAM / GREY (dividers incl)	\$ 2,879.00	\$ 1,439.50	2
WAT905370	Anaesth Cart TRIBUTE 90 - 8H DR121 w/ Digi-Lock & RHS Handle - Tower 4 & 5 Tilt Bins - GREEN (divs i	\$ 3,339.00	\$ 1,669.50	2
WAT905400	Anaesth Cart TRIBUTE 90 - 8H DR121 w/ BA Lock & RHS Handle - Tower 4 & 5 Tilt Bins - GREEN (divs inc	\$ 3,159.00	\$ 1,579.50	1
WAT907150	Cart Accessory - Aidacare/Tribute - Tilt Bin Organiser 5 Compartment	\$ 145.00	\$ 43.50	9
WAT907155	Cart Accessory - Aidacare/Tribute - Tilt Bin Organiser 6 Compartment	\$ 132.00	\$ 39.60	6
WAT907270	Avalo Handle & End Panel- 8H (K) - Light Cream/Grey	\$ 425.00	\$ 127.50	1
WAT907322	Avalo 3in Drawer (Complete) - with Light Cream/Grey Tab (K)	\$ 399.00	\$ 119.70	17
WAT907323	Avalo 6in Drawer (Complete) - with Light Cream/Grey Tab (K)	\$ 399.00	\$ 119.70	31
WAT907324	Avalo 9in (10in) Drawer (Complete) - with Light Cream/Grey Tab (K)	\$ 422.00	\$ 126.60	17
WAT907325	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Grey	\$ 2.31	\$ 0.69	111
WAT907326	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Blue	\$ 2.31	\$ 0.69	106
WAT907327	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Yellow	\$ 2.31	\$ 0.69	107
WAT907328	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Green	\$ 2.31	\$ 0.69	151
WAT907340	Drawer Sub Divider - 7in (all Drawers) - TRIBUTE	\$ 6.75	\$ 2.03	284
WAT907380	Cart Accessory - Aidacare/Tribute or AVALO - Label - Anesthesia Cart - Green	\$ 2.31	\$ 0.69	103
WAT907410	Cart Accessory - Aidacare/Tribute or AVALO - Label - Isolation Cart - Yellow	\$ 2.31	\$ 0.69	105
WAT907420	Cart Accessory - Aidacare/Tribute or AVALO - Label - Medical/Surgical Cart - Blue	\$ 2.31	\$ 0.69	103
WAT907460	Avalo Drawer Slide Assembly (K)	\$ 107.00	\$ 32.10	4
WAT907500	Avalo Keyless Lock Lever Assembly (K)	\$ 179.00	\$ 53.70	2
WAT907600	Avalo Lock Bezel Override Cover	\$ 13.80	\$ 4.14	2
WAT907620	Avalo Lock Bezel Exterior Filler - Light Cream/Grey	\$ 43.90	\$ 13.17	1
WAT907630	Avalo Aw Lock Bezel Interior Filler - Light Cream/Grey (only colour)	\$ 27.60	\$ 8.28	4
WAT907665	Avalo Writing Surface - Light Cream	\$ 173.00	\$ 51.90	1
WAT907680	Avalo Bumper - Specify Colour - Light Cream/Grey	\$ 98.00	\$ 29.40	1
WAT907682	Avalo Bumper - Specify Colour - Ext Blue	\$ 98.00	\$ 29.40	4
WAT907684	Avalo Bumper - Specify Colour - Ext Green	\$ 98.00	\$ 29.40	2
WAT907715	Avalo Keyless Override Cable	\$ 119.00	\$ 35.70	3
WAT907800	Cabinet - Imprest - 1200 x 1200 x 2000	\$ 13,629.00	\$ 4,088.70	1
WAT909230	Cart Accessory - Aidacare/Tribute - Sharp Container	\$ 184.00	\$ 55.20	2