




# GABRIEL X FIVE ACTION LIFT RECLINE CHAIR

The Aspire Gabriel X Five Action Lift Recline Chair is designed to support comfort, mobility, and ease of use. With powered adjustments for tilt, headrest, lumbar, backrest, and leg elevation, it helps users find and maintain a supported sitting position. Zoned massage and seat heating functions assist with relaxation and pain management.



- **Tilt-in-space:** Zero gravity positioning reduces forces on the lower back and improves circulation
- **Four Function Massage:** Zoned air cells alternate to relieve tight, painful lower back muscles
- **All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Electric Headrest Adjustment:** Offers customisable neck support at the press of a button
- **Adjustable Seat Heating:** Offers comfort and relief for sore muscles and joints
- **USB A and C Charging:** Located conveniently in the armrest so users can charge devices while browsing
- **Large Button Handset:** Offers high visibility back-lit buttons that are easy to use
- **Integrated Leg Extension:** Helps to comfortably support lower legs for a wide range of user heights



Zero Gravity  
Tilt-in-space



Back  
Recline



Leg  
Raise



Power  
Lumbar



Power  
Headrest



Seat  
Heating



4 Function  
Massage

DARK GREY - SMALL

EBONY - LARGE

2 YEAR LIMITED  
WARRANTY



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width
CHP227591	Dark Grey	Small	475mm	450mm	490mm	710mm	1100mm	850mm
CHP227596	Ebony	Large	500mm	520mm	530mm	780mm	1170mm	900mm