



Gabriel X5

Five Action Lift Recline Chair

Enjoy the functional finesse of our premium five-motor lift recline chair. Configure your comfort with independent adjustment of the back, leg, and headrest, combined with massage and heating for added luxury.



Zero Gravity
Tilt-in-space



Back
Recline



Leg
Raise



Power
Lumbar



Power
Headrest



Seat & Backrest
Heating



3 Function
Massage

1300 133 120 | aidacare.com.au



Zero Gravity and Tilt-in-space: Create a feeling of weightlessness by elevating your legs above your upper body, with only a limited amount of space needed behind the lift recline chair



Three Level Backrest Massage:
A choice of massage settings for personal comfort



Adjustable Seat Heating:
Integrated heating provides gentle warmth and improves comfort, perfect for year-round relaxation



Integrated Leg Extension:
Provides improved support under the lower legs for a wide range of users



Electric Headrest Adjustment:
Offers customisable neck positioning at the press of a button



Large Button Handset:
High-visibility back-lit buttons are easy to use



Stay Connected:
Conveniently placed USB-A and USB-C charging ports in both the armrest and handset for phones, tablets, and other devices



Durable and Safe:
Scratch-resistant materials and protective covers extend the life of the chair, and a built-in safety barrier keeps pets away from moving parts




Battery Backup:
Complete a lift and recline during power outages



Storage Pockets: Keep your belongings close by



Handset Loops:
Ensures the handset is always within reach



Dark Grey Microsuede



Code	Size	Seat Height	Seat Depth	Seat Width	Armrest Height	Backrest Height	Overall Height	Overall Width	Distance From Wall
CHP227592	Petite	475mm	450mm	490mm	170mm	710mm	1100mm	850mm	500mm
CHP227597	Large	500mm	520mm	530mm	170mm	780mm	1170mm	900mm	550mm