

SIGNATURE 2 DUAL ACTION LIFT RECLINE CHAIR

The Signature 2 is a dual action lift recline chair offering independent movement of the backrest and legrest for optimal positioning. Available in three size options, the Signature 2 maximises comfort and support for a wide range of users.



Mink Warwick
Global Microsuede



Lagoon Warwick
Global Microsuede



Heritage Vinyl



Custom Upholstery
Available



Adjustable Neck Pillow

CHP197781 | Mink Warwick
CHP197783 | Lagoon Warwick
CHP197782 | Heritage



Leg Raise



Back Recline

- **Three sizes help match the chair's fit to the user's body shape**, improving comfort, posture and safety when standing
- **Dual action design creates infinite positioning** via independent backrest and legrest adjustment
- **Adjustable seat depth** ensures correct positioning and safer transfers for a wide range of users
- **Integrated leg extension** helps to comfortably support lower legs for a wide range of user heights
- **Swivel castor wheels** enable easy movement of the chair for cleaning or repositioning
- **All-day battery backup** ensures continued operation during power outages, supporting safety and independence
- **Backlit hand control and integrated USB charging** are simple and convenient

2 YEAR
WARRANTY

SWL
155kg

Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228200 CHP228210 CHP228220	Mink Warwick Lagoon Warwick Heritage	Small	480mm	460 – 385mm	480mm	580mm	980mm	740mm	155kg
CHP228230 CHP228240 CHP228250	Mink Warwick Lagoon Warwick Heritage	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	155kg
CHP228260 CHP228270 CHP228280	Mink Warwick Lagoon Warwick Heritage	Large	500mm	510 – 535mm	540mm	650mm	1110mm	820mm	155kg