

Aspire Signature 2 Dual Action Lift Recline Chair

ORDER/QUOTE FORM

Submitted By:

Client/Group Name:


NDIS Participant No.

Date:

Site ID:

Sales Order No.

1. CHOOSE YOUR SIZE:




SMALL

SWL 130kg
CHP228220

\$

Base Price




MEDIUM

SWL 130kg
CHP228250

\$

Base Price







LARGE

SWL 130kg
CHP228280

\$

Base Price

2. CHOOSE YOUR FABRIC TYPE:

Swatch	Material Type	Range Name	Identify Colour
	Healthcare Vinyl	Lustrell Ranger	
	Micro Suede	Eastwood	
	Soft Touch Fabric	Manisa	
	Leather	La Casa	

3. ORDERING AND PRICING:

CHP

NEW CODE (Code + Swatch) e.g. CHP228220LT



Customer's Signature - Size, Fabric and Colour Confirmation

\$

TOTAL PRICE GST free

Confirmation Date

Once complete, please email customfurniture@aidacare.com.au

TERMS AND CONDITIONS

Custom Aspire Signature 2 Lift Recliners are made to order and therefore non-refundable. Refunds and/or exchanges will **not** be given for 'change of mind'.

By signing this form, you and/or a representative acting on behalf of the user, acknowledges that the size, fabric range and colour selection is correct.

For users with incontinence, pressure care and/or other clinical needs, it is essential that you consult your clinician, doctor or healthcare professional. This should be assessed prior to payment, to ensure that this product is right for the intended user.

For our extended list of terms and conditions, please refer to the user manual.

Processing and Payment

Full payment and/or funding application approval is required to process the order.

Lead Times

The "*Estimated Total Lead Time*" will **only** commence once full payment and/or funding approvals have been received.

*To preserve the appearance and integrity of the leather's finish, leather recliners are not made with armrest or headrest covers.