



# Headboard & Footboard Range

Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs—improving comfort for users and making life easier for carers.



Freestanding Headboards



Fixed Headboards



Footboards

## Freestanding Headboards Allure / Chesterfield

Freestanding headboards stand independently behind the bed and can be moved or removed as needed. Contemporary in design, the headboard remains stationary while the bed is raised or lowered, not impacting any wall decoration above the headboard.

- Not attached to the bed—offers full positioning flexibility
- Quick to install or remove without tools
- Softened surface that supports pillows and protects the head
- Enhances room style without limiting bed function
- Supports high-needs equipment integration in care environments
- Available in two modern styles, both at 1140mm high and 110mm deep



### Allure - Width (mm)

- [BEA830120](#) Long Single - 920
- [BEA830130](#) King Single - 1070
- [BEA830132](#) Long Double - 1370
- [BEA830134](#) Queen - 1520
- [BEA830136](#) King - 1840



### Chesterfield - Width (mm)

- [BEA830100](#) Long Single - 920
- [BEA830110](#) King Single - 1070
- [BEA830112](#) Long Double - 1370
- [BEA830114](#) Queen - 1520
- [BEA830116](#) King - 1840

## Fixed Headboards Allure

Fixed headboards attach directly to the bed frame thereby moving with the bed during tilt and height adjustments.

- By moving with the bed this reduces pillow displacement and supports the pillow and head position during tilt
- Pivoting design enhances wall protection by pivoting upon wall contact
- The Bed/Headboard can be positioned in direct proximity to a wall and still perform head-down tilt function
- Built-in handset loops (Double/Queen) keep controls within convenient reach
- Widths to match bed sizes, all 680mm high and 120mm deep

### Fixed Headboards Allure - Width (mm)

- [BEA830150](#) Long Single - 920
- [BEA830152](#) King Single - 1070
- [BEA830154](#) Long Double - 1370
- [BEA830156](#) Queen - 1520



## Footboards

The footboard completes the bed in both style and function.

- Helps complete the bed end and keep bedding in place during profiling, reducing the need for carers to continually re-adjust bedding
- Provides a practical accessory mounting surface (e.g. mattress pump)
- Portobello grey upholstery to match bed and headboards
- Secure and sturdy bolt-on attachment to bed frame
- Widths to match bed sizes, all 610mm high and 60mm deep

### Footboards - Width (mm)

- [BEA830160](#) Long Single - 920
- [BEA830162](#) King Single - 1070
- [BEA830164](#) Long Double - 1370
- [BEA830166](#) Queen - 1520

