



Montana X Maxi Single Action

1 Motor Lift Recline Chair

Stylish Bariatric Seating for Your Home

Created for comfort and bariatric support, the Montana X Maxi features a zip-access adjustable backrest for personalised postural adjustments. Its sleek design fits effortlessly into modern homes, blending clinical function with elegant style.



Effortless Control
Simple 2-button handset



**Vertical
Lift**



**Full
Recline**



**Golden Grain
Microsuede**



Trusted Bariatric Design: 220kg high safe working load and generous seating surface

Refined, Contemporary Design: Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environments

Integrated Footrest Extension: Reduces fatigue and lower limb swelling by improving healthy circulation

Protective Covers for High-wear Areas: Removable headrest and armrest covers keep your chair clean and long-lasting

Commercial-grade Seat Cushioning: Enhances comfort and durability, providing lasting support

Promotes Independent Daily Living: Simple 2-button handset assists users with cognitive or visual challenges

All-day Battery Backup: Reliable Lithium-ion battery ensures continuous operation during power outages, maintaining your safety and independence



**2 YEAR
WARRANTY**



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Depth	Overall Width
CHP227401	Golden Grain	Maxi	510mm	530mm	640mm	670mm	1060mm	940mm	1010mm