

PRODUCT UPDATE



Signature 2 Dual Action Lift Chair - Upgrades

Our Signature 2 Recliners have undergone some upgrades on the Dual Action Models.

- Integrated Leg Extension: Helps to comfortably support lower legs for a wide range of user heights. *Previously, this was a separate support piece that extended out from the legrest. The new legrest tracks the natural movement of the knee joint offering improved support when raising and lowering.*
- Custom Upholstery: Like our Raphael Chairs, Signature 2 Chairs are now able to be customised locally in Leather, Vinyl, Micro Suede and Soft Touch fabrics from Warwick.



Operational Considerations

- Product codes remain the same for our standard colours and sizes. New codes have been created for Custom models. All codes are listed for easy reference on the next page of this PAB.
- Dual Action models with the new legrest are identifiable in their cartons by the round yellow sticker as shown below. These are being added to help with stock identification during the transition process as we run down the previous Dual Action models.
- The Custom upholstery process follows the same logic as that for our Raphael chairs. The new order form is provided with this PAB.

PRODUCT UPDATES

Dual Action Signature 2 Chairs - Standard

Small	
CHP228200	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Mink - 155kg ^
CHP228210	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Lagoon - 155kg ^
CHP228220	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Heritage Vinyl - 155kg ^
Medium	
CHP228230	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Mink - 155kg ^
CHP228240	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Lagoon - 155kg ^
CHP228250	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Heritage Vinyl - 155kg ^
Large	
CHP228260	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Mink - 155kg ^
CHP228270	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Lagoon - 155kg ^
CHP228280	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Heritage Vinyl - 155kg ^

Dual Action Signature 2 Chairs – Custom

Small	
CHP228220FA	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Soft Touch Fabric (choose colour)
CHP228220LT	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Leather (choose colour)
CHP228220MS	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Microsuede (choose colour)
CHP228220VL	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Small - Healthcare Vinyl (choose colour)
Medium	
CHP228250FA	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Soft Touch Fabric (choose colour)
CHP228250LT	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Leather (choose colour)
CHP228250MS	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Microsuede (choose colour)
CHP228250VL	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Medium - Healthcare Vinyl (choose colour)
Large	
CHP228280FA	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Soft Touch Fabric (choose colour)
CHP228280LT	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Leather (choose colour)
CHP228280MS	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Microsuede (choose colour)
CHP228280VL	Chair Powerlift Recline - Aspire Signature 2 - Dual Action - Large - Healthcare Vinyl (choose colour)