

freedom **fhc**
HEALTHCARE

PRODUCT CATALOGUE

MARCH 2021



PRODUCTS THAT EMPOWER YOUR LIFE

Our customers are at the heart of everything we do. As such, our exclusive range of products have been carefully designed to empower your life by maximising your independence and quality of life.



TABLE OF CONTENTS

CHAIRS	4
LOW BACK DAY CHAIRS	HIGH BACK DAY CHAIRS
BEDROOM AIDS	12
BEDSIDE COMMODES	BED STICKS AND POLES
BATHROOM AIDS	18
SHOWER STOOLS AND CHAIRS	BATH TRANSFER BENCH
OVER TOILET AIDS AND SURROUNDS	TOILET SEAT RAISER
MOBILITY AIDS	34
ROLLATORS	WALKING FRAMES
FOREARM WALKERS	CRUTCHES AND WALKING STICKS




CHAIRS

LOW BACK DAY CHAIRS

HIGH BACK DAY CHAIRS












LSR535 LOW BACK DAY CHAIR WITH OVAL TUBING




-  BURGUNDY VINYL
CHP208300
-  BLUE VINYL
CHP208310
-  RUSTIC VINYL
CHP208320







-  ARTG 140476
-  AS/NZS 4688.3:2000
-  AS/NZS 4688.2:2000

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	UNIT WEIGHT	WEIGHT CAPACITY
 880-1030mm	 585mm	 620mm	 480mm	 485mm	 465-615mm	 185mm	 8kg	 180kg

BC1 LOWBACK DAY CHAIR FABRIC AND VINYL



-  DENIM FABRIC
CHP196700
-  LATTE VINYL
CHP196701

SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	PRODUCT WEIGHT	SWL
 490mm	 535mm	 380-530mm	 550mm	 15kg	 160kg

BC1 HIGH BACK DAY CHAIR

FABRIC AND VINYL



DENIM FABRIC
CHP196500

LATTE VINYL
CHP196501

SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	PRODUCT WEIGHT	SWL
470mm	535mm	380-530mm	630mm	16kg	160kg

BC2 WATERFALL DAY CHAIR

FABRIC AND VINYL




CHARCOAL FABRIC
CHP196799

INK VINYL
CHP196800

SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	PRODUCT WEIGHT	SWL
780mm	525mm	430-580mm	650mm	16.5kg	160kg

BC4 STANDARD

FABRIC AND VINYL



INK VINYL
CHP196810

LATTE VINYL
CHP196811

PRODUCT CODE: CHP196810

SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	PRODUCT WEIGHT	SWL
440-510mm	510mm	440-590mm		16kg	160kg



**SHOP THE FULL
FREEDOM HEALTHCARE
RANGE ONLINE TODAY!**

[Home](#) | [Bathroom](#) | [Bedroom](#) | [Mobility & Access](#) | [Seating](#)



TO SCAN QR CODE

- ① Open phone camera
- ② Hover over QR code
- ③ Click on link at the top




BEDROOM AIDS

BEDSIDE COMMODES

BED STICKS & POLES



JBA500-C BEDSIDE COMMODE WITH LUXURY PADDED UPHOLSTERY



BEIGE

PRODUCT CODE: BEC046405

ARTG 135307

AS/NZS 3973:2009

12 MONTH WARRANTY

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	UNIT WEIGHT	WEIGHT CAPACITY
935-1085mm	600mm	620mm	490mm	48.5kg	500-650mm	170mm	10.1kg	160kg

JBA506 HEIGHT ADJUSTABLE BED GRAB BAR ALUMINIUM



Designed to assist the patient into a seated position

Can be used on either side of the bed



PRODUCT CODE: BEA831100


OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	UNIT WEIGHT
420-550mm	550mm	830mm	1.9kg

ARTG 174584

AS EN 12182:2015

12 MONTH WARRANTY


JBA503 ADJUSTABLE BED CRADLE TO KEEP BEDDING OFF THE PATIENT



Designed to keep the bedding elevated off the patient





Height adjustable

Sturdy powder coated steel frame



Frame fits securely under the mattress


PRODUCT CODE: BEA831110

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	UNIT WEIGHT
 450-675mm	 710mm	 360mm	 2.8kg

12 MONTH WARRANTY

✓ ARTG 174584 ✓ AS EN 12182:2015

JBA505-KD COBRA BED STICK SUPPORT POLE



Designed to assist the patient into a seated position





Powder-coated aluminium frame

Base fits securely under mattress



Can be used on either side of the bed

PRODUCT CODE: BEA831120

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	UNIT WEIGHT
 500mm	 390mm	 470mm	 1.3kg

12 MONTH WARRANTY

✓ ARTG 174584 ✓ AS EN 12182:2015



BATHROOM AIDS

SHOWER STOOLS & CHAIRS

OVER TOILET AIDS & SURROUNDS

BATH TRANSFER BENCH

TOILET SEAT RAISER



HBA370 SPLIT SEAT CHAIR - 3 in 1

OVER TOILET AID - SHOWER CHAIR - COMMODE



HBA406 OVER TOILET AID

OVAL TUBE ALUMINIUM



OPTIONAL EXTRAS




BACKREST HBA371
(BTT163120)



SPLASH GUARD TSA370-B
(GSP650720)



EXTENSION LEGS TSA370-K
(GSP650100)



OPTIONAL EXTRAS



EXTENSION LEGS TSA370-K
(GSP650100)



ARMREST PAD TSA370-M
(GSP650460)



ASPIRE COMMODE BOWL WITH HANDLE & LID
(BTC081610)

12 MONTH WARRANTY

PRODUCT CODE: BTT163100

ARTG 135307

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
650-800mm	595mm	645mm	450mm	480mm	440-590mm	210mm	6.3kg	160kg

12 MONTH WARRANTY

PRODUCT CODE: BTT163140

ARTG 135307

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
650-800mm	605mm	615mm	430mm	460mm	440-590mm	210mm	3.8kg	200kg

HBA407 SHOWER CHAIR

OVAL TUBE ALUMINIUM




OPTIONAL EXTRAS



EXTENSION LEGS TSA370-K
(GSP650100)



ARMREST PAD TSA370-M
(GSP650460)



HBA408 SHOWER STOOL

OVAL TUBE ALUMINIUM



OPTIONAL EXTRAS



EXTENSION LEGS TSA370-K
(GSP650100)



BACKREST TSA408-C
(BTS150125)



ARMREST PAD TSA370-M
(GSP650460)

12 MONTH WARRANTY

ARTG 135308

AS/NZS 3973:2009

PRODUCT CODE: BTS150130

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
795-945mm	605mm	615mm	400mm	460mm	440-590mm	210mm	3.8kg	200kg

12 MONTH WARRANTY


ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
670-820mm	520mm	610mm	290mm	480mm	480-630mm	190mm	4.2kg	200kg

HBA424 BATH TRANSFER BENCH

HEAVY DUTY OVAL TUBE ALUMINIUM



HBA402 SHOWER CHAIR - Premium


ALUMINIUM



OPTIONAL EXTRAS



EXTENSION LEGS TSA370-K (GSP650100)



12 MONTH WARRANTY

PRODUCT CODE: BTS153100

ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	TUB CLEARANCE HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
825-975mm	520mm	785mm	410mm	680mm	490-640mm	275mm	6kg	200kg

12 MONTH WARRANTY

PRODUCT CODE: BTS150110

ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
790-940mm	640mm	600mm	410mm	490mm	440-590mm	210mm	3.3kg	140kg

HBA391-A SHOWER CHAIR ALUMINIUM



12 MONTH WARRANTY

PRODUCT CODE: BTS150200

ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
770-895mm	510mm	550mm	400mm	470mm	415-540mm	190mm	3.6kg	130kg

HBA403 SHOWER STOOL - Premium ALUMINIUM



OPTIONAL EXTRA



BACKREST TSA403-C (BTS150125)

12 MONTH WARRANTY

PRODUCT CODE: BTS150120

ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
670-820mm	440mm	580mm	220mm	500mm	490-640mm	210mm	2.8kg	140kg

HBA392-A SHOWER STOOL ALUMINIUM



12 MONTH WARRANTY


PRODUCT CODE: BTS150210

ARTG 135308

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
670-795cm	390mm	490mm	300mm	495mm	460-585mm	210mm	3.4kg	130kg

HBA400 OVER TOILET AID - Premium ALUMINIUM



OPTIONAL EXTRA



ASPIRE COMMODE BOWL WITH HANDLE & LID
(BTC081610)

12 MONTH WARRANTY

PRODUCT CODE: BTT163130

ARTG 135307

AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
650-800mm	610mm	580mm	430mm	500mm	440-590mm	210mm	3.65kg	140kg

HBA390A OVER TOILET AID ALUMINIUM



HBA418 TOILET PAN SURROUND ALUMINIUM



12 MONTH WARRANTY

PRODUCT CODE: BTT163200

✓ ARTG 135307

✓ AS/NZS 3973:2009

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
605-730mm	510mm	550mm	370mm	470mm	410-530mm	190mm	4kg	130kg

PRODUCT CODE: BTT165110

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
660-785mm	470mm	540mm	2kg	120kg

12 MONTH WARRANTY

✓ ARTG 135307

✓ AS/NZS 3973:2009

HBA425 BATH TRANSFER BENCH WITH ROTATING SEAT



12 MONTH WARRANTY

ARTG 135308 AS/NZS 3973:2009

PRODUCT CODE: BTS153110

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	TUB CLEARANCE HEIGHT	ARMREST HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
930-1030mm	510mm	800mm	370mm	430mm	400-500mm	100mm	8.8kg	130kg

HBA444 TOILET SEAT RAISER WITH FLIP BACK ARMREST



PRODUCT CODE: BTT164100

OVERALL HEIGHT	OVERALL DEPTH	OVERALL WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
50, 100 or 150mm	430mm	380mm	3.8kg	136kg

12 MONTH WARRANTY
ARTG 135307 AS/NZS 3973:2009

MOBILITY AIDS

WHEELCHAIRS

ROLLATORS

FOREARM WALKERS

WALKING FRAMES

CRUTCHES & WALKING STICKS



AWC132-133 EXCEL SUPERLITE SELF-PROPELLED WHEELCHAIR



Spillproof nylon fabric upholstery

220mm Armrest Height

Double-braced aluminium construction

Quick release removable spoke wheels

Airless "soft ride" polyurethane tyres

Folding backrest for easy transportation

Flip Up Removable Padded Armrest

Park brake below seat level when engaged for easy side transfer

Outwards Swinging Removable Legrest

300-400mm Footrest Length

Heel Straps

150mm Solid Front Castors

Height Adjustable Plastic Footplates

OPTIONAL EXTRA



300mm TRANSIT WHEELS
TSA117AE

BLUE
MWS452180 - AWC132 (510mm)
MWS452170 - AWC133 (460mm)

12 MONTH WARRANTY

✓ ARTG 141450 ✓ ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	FOLDED WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
950mm	1070mm	AWC132: 690mm AWC133: 640mm	430mm	AWC132: 510mm AWC133: 460mm	500mm	450mm	290mm	12.7kg	110kg

AWC136-137 EXCEL SUPERLITE TRANSPORTER TRANSIT WHEELCHAIR



Handbrake with park locking mechanism connected to hub brakes

Folding backrest for easy transportation

Double-braced aluminium construction

300mm Rear Wheels

Drum Hub Brake System

Airless "soft ride" polyurethane tyres on low maintenance mag wheels

Flip Up Removable Padded Armrest

Spillproof Nylon Fabric Upholstry

Outwards Swinging Removable Legrest

300-400mm Footrest Length

Heel Straps

150mm Solid Front Castors

Height Adjustable Plastic Footplates

BLUE
MWS452120 - AWC136 (510mm)
MWS452110 - AWC137 (460mm)

12 MONTH WARRANTY

✓ ARTG 141450 ✓ ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	FOLDED WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
950mm	920mm	AWC136: 690mm AWC137: 640mm	430mm	AWC136: 510mm AWC137: 460mm	500mm	450mm	290mm	11.8kg	110kg

BRO204-HD GENESIS 3 - STANDARD HEAVY DUTY ROLLATOR



BRO209-HD GENESIS 3 - EXTRA WIDE HEAVY DUTY ROLLATOR



Rear Frame widened for maximum walking space



Larger wider wheels minimise wobble and increase grip.



Oval Tube Construction for style and strength



Rear Frame widened for maximum walking space



Larger wider wheels minimise wobble and increase grip



Oval Tube Construction for style and strength

RED
WAF750310

CHAMPAGNE
WAF750280

12 MONTH WARRANTY

ARTG 141450 ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
900-1020mm	770mm	680mm	330mm	490mm	550mm	380mm	8 inch	9.4kg	270kg

BLUE
WAF750335

12 MONTH WARRANTY

ARTG 141450 ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
900-1020mm	770mm	750mm	360mm	550mm	550mm	380mm	8 inch	9.9kg	270kg

BRO201-028 8" QUAD WALKER

ROLLATOR WITH ADJUSTABLE SEAT HEIGHT



FOLDING WALKER
For easy transport and storage

BRO202-02 MODULA - BASIC

ROLLATOR



RED
WAF750530

ARTG 141450

ISO 11199-2:2005

12 MONTH WARRANTY

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
770-960mm	695-780mm	620mm	330mm	440mm	500-580mm	320mm	8 inch	7.8kg	140kg

RED
WAF750430

BLUE
WAF750380

PURPLE
WAF750410

ARTG 141450

ISO 11199-2:2005

12 MONTH WARRANTY

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
820-970mm	690mm	620mm	330mm	440mm	560mm	230mm	8 inch	8.4kg	130kg

BRO202-02-Y MODULA - DELUXE ROLLATOR













Rear Frame widened for maximum walking space

- RED WAF750470
- BLUE WAF750440
- PURPLE WAF750441

12 MONTH WARRANTY

- ARTG 141450
- ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
									
820-940mm	690mm	620mm	330mm	440mm	550mm	230mm	8 inch	8.4kg	140kg

BRO199 TRANSROLLER COMBO (TRANSPORTER/ROLLATOR)



Double-X Frame for additional strength













Swing away/removable legrests

Height adjustable footrests

- SILVER WAF750340

12 MONTH WARRANTY

- ARTG 141450
- ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
									
860-980mm	roll: 710mm stroll: 980mm	640mm	310mm	460mm	560mm	270mm	8 inch	9.3kg	100kg

BRO200S MOBYFLEX FOLDING ROLLATOR



MOBYFLEX FOLDING SYSTEM



FOLDING REAR LEGS



FOLDING FRONT WHEELS



FULLY FOLDED
Ultra-compact for easy transportation.

ATTACHMENTS



WALKING STICK HOLDER
Attachment Included

STRIDE WALKER ROLLATOR



12 MONTH WARRANTY

ISO 11199-2:2005

SILVER WAF750360

ARTG 141450

ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
780-900mm	690mm	640mm	250mm	460mm	530mm	270mm	8 inch	7.6kg	130kg

RED WAF750150

ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
790-980mm	830mm	600mm	230mm	480mm	520mm	230mm	Front - 11 inch Back - 9 inch	9kg	135kg

BRO212A-HD FOREARM WALKER

HEAVY DUTY



BRO212B-HD FOREARM WALKER

HEAVY DUTY



SILVER
WAF699120

12 MONTH WARRANTY
ARTG 141450 ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	ARMREST LENGTH	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
1135-1355mm	860mm	690mm	430-530mm	260mm	8 inch	10.6kg	220kg

SILVER
WAF699201

12 MONTH WARRANTY
ARTG 141450 ISO 11199-2:2005

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	ARMREST LENGTH	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
1135-1355mm	880mm	660mm	360-460mm	270mm	Front: 8 inch Rear: 5 inch	9.3kg	220kg

BRO213-BHD FOREARM WALKER WITH SEAT



BRO214 VIKING FOREARM WALKER

DESIGNED FOR HOSPITALS AND OTHER HEALTH CARE INSTITUTIONS



BLUE WAF699210		12 MONTH WARRANTY								
OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	ARMREST LENGTH	FOLDED WIDTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
930-1160mm	770mm	680mm	330mm	490mm	550mm	320-410mm	380mm	8 inch	9.2kg	200kg

BLUE WAF699130		12 MONTH WARRANTY				
OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	CASTOR HEIGHT	WEIGHT CAPACITY	ARTG 141450	ISO 11199-2:2005
920-1350mm	780mm	715-915mm	4 inch	200kg	ARTG 141450	ISO 11199-2:2005

BRO215 TITAN FOREARM WALKER

DESIGNED FOR HOSPITALS AND OTHER HEALTH CARE INSTITUTIONS



CWF230-A QUICK FOLD WALKING FRAME

ALUMINIUM



OPTIONAL EXTRAS

ACCESSORIES PACK
TSA215 - H
WAF699216



OPTIONAL EXTRAS



Folds compact for storage and transportation.

12 MONTH WARRANTY

ARTG 141450 ISO 11199-2:2005

BLUE
WAF699215

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	ARMREST LENGTH	CASTOR HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
950-1280mm	840mm	705mm	380-510mm	4 inch	8kg	200kg

PRODUCT CODE: WAF752100

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	OPTIONAL CASTOR	OVERALL WEIGHT	WEIGHT CAPACITY
790-1000mm	640-780mm	570mm	210mm	200mm	6.9kg	130kg

12 MONTH WARRANTY

ARTG 141450 ISO 11199-2:2005

CWF233 ULTRALITE WALKING FRAME

ALUMINIUM



Folds compact for easy storage and transportation.



PRODUCT CODE: WAF752140

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
810-910mm	470mm	580mm	100mm	2.3kg	100kg

ARTG 141450

12 MONTH WARRANTY

ISO 11199-2:2005

CWF234 DUAL HANDLE WALKING FRAME

ALUMINIUM



OPTIONAL EXTRAS



FRONT WHEELS
TSA230-F
(WAF752170)



REAR SKI TIPS
TSA230-G
(WAF752150)



PRODUCT CODE: WAF752090

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
830-1010mm	520mm	580mm	95mm	2.9kg	130kg

ARTG 141450

12 MONTH WARRANTY

ISO 11199-2:2005

RBE804 BARIATRIC WALKING FRAME

STEEL



PRODUCT CODE: WAF752162

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	OVERALL WEIGHT	WEIGHT CAPACITY
760-1010mm	460mm	740mm	95mm	3kg	270kg


✓ ARTG 141450

✓ ISO 11199-2:2005

12 MONTH WARRANTY


CWF215 UNDERARM CRUTCHES

ADJUSTABLE PLATINUM



ADJUSTABILITY RANGE

The Platinum crutch can be adjusted to match the heights of traditional crutch sizes, making it the all rounder crutch!



PRODUCT CODE: WAC692505

ARMPIT PAD HEIGHT	ARMPIT TO HANDLE	UNIT WEIGHT	WEIGHT CAPACITY
950-1540mm	230-350mm	0.9kg	120kg

✓ ARTG 227868

✓ ISO 11199-2:2005

12 MONTH WARRANTY

DQC245-247 QUAD CANES ALUMINIUM

ALL QUAD CANES FEATURE:

- Ergonomic handles made from durable plastic for a long lifespan.
- Swivel Handles to change to left or right handed
- Strengthened light-weight aluminium frame
- Height adjustable extension
- **REPLACEABLE RUBBER TIPS**
TSA245-D - 12mm (WAS748340)
TSA245-D - 14mm (WAS748350)



SMALL FLAT BASE CANE
DQC245
(WAS748321)
with small base.



PYRAMID QUAD CANE
DQC246
(WAS748320)
with larger base for greater stability.



LARGE FLAT BASE CANE
DQC247
(WAS748310)
with wider base for greater stability.

12 MONTH WARRANTY

✓ ARTG 141452 ✓ ISO 11199-1:1999

MODEL	OVERALL HEIGHT	BASE LENGTH	BASE WIDTH	BASE HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
Small Flat Base Canes (DQC245)	665-895mm	210mm	155mm	72mm	0.86kg	110kg
Pyramid Quad Canes (DQC246)	760-960mm	315mm	210mm	265mm	0.8kg	
Large Flat Base Canes (DQC247)	770-980mm	310mm	225mm	75mm	1.05kg	

EWS258-269 WALKING CANES ALUMINIUM



12 MONTH WARRANTY

✓ ARTG 141452

✓ ISO 11199-1:1999

OVERALL HEIGHT	WEIGHT CAPACITY
770-970mm	110kg





CHOOSE YOUR STYLE	SWAN NECK EWS258	CROOK HANDLE EWS259	T HANDLE PLASTIC EWS268	T HANDLE ERGO PLASTIC EWS269
BRONZE	WAS748370	WAS748160	WAS748540	WAS748490
BLACK	WAS748360	WAS748140	WAS748150	WAS748480
CHROME	WAS748380	WAS748390	WAS748550	WAS748500
DARK FLORAL	WAS748400	WAS748180	WAS748560	WAS748510
TARTAN	WAS748420	WAS748210	WAS748580	WAS748530
BEIGE LEAVES	WAS748410	WAS748190	WAS748570	WAS748520

EWS256 FOLDABLE WALKING CANE ALUMINIUM




-  **BRONZE**
EWS256-01
(WAS748240)
-  **BLACK**
EWS256-02
(WAS748230)







OVERALL HEIGHT	WEIGHT CAPACITY
 800-900mm	 114kg

12 MONTH WARRANTY
 ARTG 141452  ISO 11199-1:1999

RBE806 BARIATRIC QUAD CANE REINFORCED STEEL




PRODUCT CODE: WAS748290

OVERALL HEIGHT	BASE LENGTH	BASE WIDTH	BASE HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
 710-930m	 240mm	 180mm	 80mm	 1.5kg	 270kg




12 MONTH WARRANTY
 ARTG 141452  ISO 11199-1:1999

RBE808 BARIATRIC WALKING CANE

REINFORCED STEEL



PRODUCT CODE: WAS748100

OVERALL HEIGHT	OVERALL WEIGHT	WEIGHT CAPACITY
 710-930mm	 0.9kg	 270kg

12 MONTH WARRANTY

✓ ARTG 141452 ✓ ISO 11199-1:1999

TO SEE MORE, VISIT OUR WEBSITE

freedomhealthcare.com.au



TO SCAN QR CODE

- ① Open phone camera
- ② Hover over QR code
- ③ Click on link at the top

COPYRIGHT

© Copyright Aidacare Pty Ltd 2021 Version 1.0

No part of this catalogue may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means.

DISCLAIMER

Due care has been taken to ensure all the information contained in this catalogue is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. We reserve the right to modify the design or appearance of any product displayed in this catalogue without prior notice.



www.freedomhealthcare.com.au