

teal

Myla Security Bed

A secure freestanding bed designed for challenging environments. Smooth, round profile edges enhances safety and security.



**Colour shown for illustrative purposes only.*

- Durable MDF vinyl wrapped soft edge top (25mm thick)
- Channeled platform allows for air to circulate under mattress
- Waterfall edging around the entire perimeter prevents fluid ingress and facilitates easier cleaning for improved infection control
- Sloping sides enable easy staff access for cleaning or assisting a user to aid sit-to-stand movements
- Integrated mattress retainers prevent mattress movement
- Tested and compliant to BS EN 597-1, BE EN 597-2, BS6807:2006 clause 9 & BS 7177:2008



Large range of mattresses available
Contact us to help find your solution



Discover a spectrum of colours | [Contact us for more information](#)

BEB050500
Myla Security Bed

Overall Height	Overall Width	Overall Depth	Suggested Mattress Size	Product Weight
385mm	1020mm	1950mm	1900 x 915mm	78kg

Exclusively distributed by
 **aidacare**
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au